

PHILIPS

USB cable

Philips Avent


Breast pumps

Compatible with only USB
charged devices

White

Cord length: 125 cm

CP2212



A USB cable to charge your device

Check for compatibility below

A USB cable to effectively charge your Hands-free pump. Only compatible with the Philips Avent Hands-free electric breast pump. Use the USB cable with a 5 V/2 A or 10 W adapter or other power supply (power banks, USB wall sockets etc.).

Specifications

Replaceable part

Fits product types: SCD553, SCF531, SCF532,
SCF533, SCF534, SCF537, SCF539, SCF543, SCF547

