

PHILIPS

AVENT



Single Electric Breast
Pump

Hands-free



SCF543/11

Free your hands. Free your time.

Mimics your baby's drinking rhythm

The Philips Avent Hands-free Single Electric Breast Pump is made for effective* and comfortable expression. It mimics your baby's rhythm and with power, gentleness and lightness, it's your ideal hands-free pump companion.

Physical freedom while expressing

- Move freely with the ultra-light collection cup
- Wearable comfort for maximum movement

Mimics baby's drinking rhythm

- Express effectively* with a rhythm that mimics your baby's
- Quiet hospital-strength pumping

Easy to use, simple to clean and assemble

- Track milk flow with the fully transparent collection cup
- Personalise your pumping session
- Easy to clean and assemble with few parts

Inclusive, seamless fit

- A range of breast shield and insert sizes for every mum
- Silicone shield gently shapes to your breast

Hands-free

Single Electric Breast Pump

SCF543/11

Highlights

Mimics baby's drinking rhythm



Babies know how to drink best! That's why our hands-free electric breast pump mimics their natural drinking rhythm – so you can find your own rhythm for effective* expression and optimal milk output. Pumping up to 85 times per minute, the Philips AVENT Hands-free Breast Pump is 2 x faster than most other pumps**.

Hospital-strength motor




The nearly silent motor unit delivers the hospital-strength power you expect from a traditional breast pump, yet from a compact and wearable unit. With its rechargeable battery, the cordless motor unit pumps powerfully and discreetly anywhere.

Unique SkinSense shields



Our unique SkinSense silicone breast shield is designed for ultimate comfort, softly shaping around your breast using your body's natural warmth. The breast shield is made with food-grade silicone for complete peace of mind.

Transparent collection



The see-through collection cup provides an unobstructed view inside. It helps you correctly place your nipple for your best milk output and visually track the volume of breast milk. For confident pumping from the first drop until finish.

Ultra-light collection cup




Our collection cup weighs just a little over 100 grams when empty - up to 3 x lighter than a fully wearable breast pump. The cup sits comfortably in-bra with barely a sign that it's there.

Wear it your way



Our pump comes with a braided carrying cord and apparel clip for wearing the motor unit cross-body or attached to your waistband. With the freedom to wear it your way, you stay in your flow knowing control is right by your side. Use the tube coupler and cord adjuster to precisely fine-tune your fit.

Personalised pumping




Get the most from every session with 8 suction levels for stimulating the flow of milk and 16 suction levels for expressing it. You can fine-tune as you pump and with your ideal setting automatically memorised, your pump is always ready to go.

Easy to clean and assemble



We designed our pump and accessories to be as simple as possible. With only a few connecting pieces, each collection cup is a breeze to take apart, clean (by hand or in a dishwasher), disinfect and put back together.

Find the right fit



We support every mum to feel their best and having the right fit is an essential part of that. Our range of breast shields and shield insert sizes fits 99% of mums. You can experience your ultimate comfort while pumping. Available to buy separately.

Hands-free

Single Electric Breast Pump

SCF543/11

Specifications

Power

Source: Rechargeable battery / USB-C (adapter not included)

Material

Cup, front cap, breast shield cover:

Polypropylene

Breast shield, membrane, valve, insert: Liquid silicone rubber

Silicone tube: Silicone

What is included

Breast milk bags: 25 pcs

Carrying cord: 1 pcs

Rechargeable motor unit (with a detachable

clip): 1 pcs

Collection cup (assembled): 1 pcs

Silicone tubes: 1 pcs

Breast shield (dust) cover: 1 pcs

Breast shield (21 mm, 24 mm): 1 pcs each

125 ml (4 oz) Natural Response bottle (with

Flow 2 teat): 1 pcs

Nipple sizing card: 1 pcs

Insert 19 mm: 1 pcs

USB cable: 1 pcs

Ease of use

Breast pump use: Easy to use, simple to clean and assemble

Functions

Flexibility in movement: Hands-free pumping

Personalisation settings: Stimulation mode with 8 suction levels and expression mode with 16 suction levels

Motor unit: Rechargeable battery / USB-C, Fully charging the battery takes up to 1.5 hours and lasts approx. 5 pumping sessions

Breast shield and insert specifics: Several sizes available for the best fit

