

PHILIPS

Silicone tube

Philips Avent

Breast pumps

Transparent

1 pc in a pack

CP2211



Connecting the collection cups with the motor unit

Check for compatibility below

The silicone tube connects the Philips Avent Hands-free electric collection cups with the motor unit. Do not trim the tubes as this increases the suction level, which may hurt or feel uncomfortable and lower your milk production.

Specifications

Replaceable part

Fits product types: SCD553, SCF439, SCF442,
SCF481, SCF531, SCF532, SCF533, SCF534,
SCF537, SCF539, SCF543, SCF547

