

# PHILIPS

Silicone diaphragm


Philips Avent

Breast pumps

Pink

2 pc in a pack

CP2218



## Creates the vacuum needed for expression

Check for compatibility below

This diaphragm creates the vacuum needed for expressing breast milk. For use with Philips Avent Hands-free electric breast pumps.

# Specifications

## Replaceable part

Fits product types: SCD553, SCF439, SCF442,  
SCF481, SCF531, SCF532, SCF533, SCF534,  
SCF537, SCF539, SCF543, SCF547

