

# PHILIPS

## sonicare



Sonic electric  
toothbrush

4100 Series

Sonicare technology

QuadPacer and SmarTimer

Slim ergonomic design

14-day battery life



HX3689/23

## Goodbye manual toothbrush. Hello Sonicare.

Up to 7X better plaque removal\*

Sonicare technology combined with our brush action gently removes up to 7x better\* plaque than a manual toothbrush.

### **Designed around you**

- Use less. Reuse more
- Always know when to replace your brush heads
- 14-day long lasting battery life
- Switch between 2 intensity settings

### **Optimize your brushing**

- Optimize your brushing with SmarTimer and QuadPacer

### **Protects your gums**

- Pressure Sensor helps protect your teeth and gums

### **Provides superior clean**

- Our unique technology brings you a powerful yet gentle clean
- Up to 7x better plaque removal vs a manual toothbrush\*


## 4100 Series

Sonic electric toothbrush

HX3689/23


# Highlights

### Removes plaque 7x better\*



Sonicare electric toothbrush with advanced sonicare technology is clinically proven to remove up to 7x better\* plaque from your teeth and along the gum line while protecting your gums.

### Protects your gums



Built-in pressure sensor automatically detects pressure you apply, warns you, and reduces the toothbrush vibrations automatically to help protect your gums. The toothbrush will make a pulsing sound as a reminder to ease off the pressure. 7 out of 10 people found this feature helped them become better brushers.

### SmartTimer and QuadPacer



2 minute SmartTimer and 30-second QuadPacer guide you to brush for the recommended time in all areas of your mouth to ensure a complete clean.

### 2 intensity settings



Switch between 2 intensity settings to make the brushing session comfortable for you.

### Advanced Sonicare Technology



Powerful bristle vibrations drive microbubbles deep between your teeth and along your gum line for a refreshing experience. You will get two month's worth of manual brushing in just 2 minutes.\*\* 31000 brush strokes per minute gently clean your teeth, break up plaque, and sweep it away for an exceptional daily clean.

### BrushSync replacement reminder



All brush heads wear out over time, so you'll want to keep an eye on yours to make sure you're still getting a great clean. Our BrushSync technology tracks how long you've been using your brush head for, and how hard you've been

brushing. The BrushSync replacement reminder on your handle and a short beep gives you a heads up when it's time to replace it.

### Long lasting battery life



Battery life lasts up to 14 days so you can go long periods without charging.

### Pulling the plug on waste



Our ambition at Philips is to reduce electronic waste. One way we do this is to minimize the number of adapters we put in the market. These products come without a power adapter because we are committed to reducing electronic waste.

## 4100 Series

Sonic electric toothbrush

HX3689/23

# Specifications

### Power

Voltage: DC5V

### Technical specifications

Battery: Rechargeable

Operating time (full to empty): 14 days\*\*\*

Battery type: Lithium ION

### Design and finishing

Color: Cleansing Blue Prismatic

### Service

Warranty: 2-year limited warranty

### Ease of use

Battery indicator: Light shows battery status

Brush head system: Easy click-on brush heads

Handle: Slim ergonomic design

BrushSync Replacement Reminder: Replacement reminder, Always know when to, replace brush heads, Reminder icon lights up, to always ensure best results

### Items included

Charger: 1

Brush heads: 1 C2 Optimal Plaque Control

Handle: 1 4100 Series

### Cleaning performance

Performance: Removes up to 7x more plaque\*

Timer: Quadpacer and SmarTimer

Pressure feedback: Vibrates handle to alert user

### Modes

Clean: For exceptional everyday clean

2 intensities: High, Low

### Smart sensor technology

BrushSync Replacement Reminder: Always know when to, replace brush heads

Pressure sensor: Alerts when brushing too hard

