

PHILIPS

sonicare

Sonic electric
toothbrush

ProtectiveClean 4300

Built-in pressure sensor

1 cleaning mode

1 x BrushSync feature



HX6807

Cleaner teeth. Gently does it.

Remove up to 7 x more plaque than a manual

Feel the difference of a gentle clean with our Pressure Sensor while removing up to 7 x more plaque than a manual toothbrush

Proven to improve oral health

- Remove up to 7x more plaque with our InterCare brush head*
- Safe and gentle on sensitive areas, orthodontics and dental work

Designed around you

- A great way to get used to your new toothbrush
- Encouragement to brush thoroughly

Innovative technology

- Always know when to replace your brush heads
- Connects smart brush handle and smart brush heads
- Philips Sonicare's advanced sonic technology
- Lets you know when you're pressing too hard

Highlights

Remove up to 7 x more plaque



Click on our InterCare brush head to improve gum health in just 2 weeks. Extra-long bristles help remove more plaque in hard-to-reach spots and between teeth for healthy gums.

Safe and gentle



You can be sure of a safe brushing experience: our sonic technology is suitable for use with braces, fillings, crowns and veneers, and helps prevent cavities and improve gum health.

Philips Sonicare technology



Philips Sonicare's advanced sonic technology pulses water between your teeth, and its brush strokes break plaque up and sweep it away for an exceptional daily clean.

Too much pressure?



You may not notice if you're brushing too hard, but pressure sensor on your electric toothbrush will. If you start taking it out on your gums, the toothbrush will make a pulsing sound as a reminder to ease off on the pressure.

BrushSync technology



A microchip-enabled technology that detects and synchronises the smart brush head with the smart handle. The smart handle and smart brush head pair is a powerful combination that enables Smart replacement reminders.


BrushSync replacement reminder



All brush heads wear out over time. But our BrushSync technology tracks how long you've been using your brush head for, and how hard

you've been brushing. When it's time to replace it, a light on your handle and a short beep will let you know. That way, you can be sure your brush head is doing a good job.

Easy-start programme



Starting something new can take getting used to. That's why our easy-start programme gives you the option of a gradual, gentle increase in brushing power the first 14 times you use your new toothbrush.

Helpful timers



Our QuadPacer lets you know when you've spent just the right amount of time cleaning each part of your mouth, while our Smartimer tells you when you've brushed for the recommended two minutes.

ProtectiveClean 4300

Sonic electric toothbrush

HX6807/63

Specifications

Power

Voltage: 110-220 V

Technical specifications

Operating time (full to empty): Up to 2 weeks

Battery: Rechargeable

Battery type: Lithium ION

Energy consumption: Standby without display
<0.5 W

Design and finishing

Colour: White and Mint

Service

Warranty: 2-year limited warranty

Ease of use

Handle compatibility: Easy click-on brush heads

Handle: Slim ergonomic design

Battery indicator: Light shows battery status

Brushing time: Up to 2 weeks

Items included

Handle: 1 ProtectiveClean

Brush heads: 1 i InterCare standard

Charger: 1

Cleaning performance

Performance: Removes up to 7x more plaque*

Health benefits: Helps reduce cavities

Timer: Quadpacer and Smartimer

Speed: Up to 62,000 brush movements/min

Pressure feedback: Vibrates handle to alert user

Modes

Clean: For exceptional everyday cleaning

2 intensities: Low, High

Smart sensor technology

Pressure sensor: Alerts when brushing too hard

BrushSync Replacement Reminder: Always know when to, replace brush heads

