



PHILIPS

Philips Easy Scan Plan

Drive CMR access by simplifying cardiac planning

Easy Scan Plan is a planning solution that significantly simplifies the cardiac planning workflow. With one 3D survey scan you can view and plan the required cardiac orientations. It is especially useful for complex procedures, such as double-oblique planning. Easy Scan Plan offers simplified cardiac planning for technologists and reduces the number of breath holds required from the patient. Additionally, it is a versatile solution that can be used for all cardiac patients, including complex cases like congenital heart disease.

Product benefits

- 1 single 3D survey scan
- View required cardiac orientations while planning
- Easy cardiac planning workflow

Field strength	1.5T, 3.0T
Main applications	Cardiac
Sequence	Cardiac survey scans

Philips Easy Scan Plan



View the required cardiac orientations during planning.

