

# PHILIPS

## AVENT



### Double Collection Cups

### Hands-free

Mimics baby's drinking rhythm

Freedom while expressing

Easy to use, clean & assemble

Inclusive, seamless fit



SCF439/01

## Convenient for home and on the go

### Collection cups for hands-free pumping

Whether at home, on the go, or even back at work, having an extra pair of hands-free collection cups ready can make your pumping routine simpler and more efficient. Includes 4x SkinSense silicone breast shields, and 2x shield inserts.

#### Other benefits

- For Premium, Advanced and Hands-Free electric breast pumps

#### Physical freedom while expressing

- Move freely with the ultra-light collection cups

#### Mimics baby's drinking rhythm

- Effective\* expression with a rhythm that mimics your baby's

#### Easy to use, simple to clean and assemble

- Track milk flow with the fully transparent collection cups
- Easy to clean and assemble with few parts

#### Inclusive, seamless fit

- The complete range of breast shields and insert sizes
- Find your perfect fit in seconds
- Silicone shields gently shape to your breasts

## Hands-free

Double Collection Cups

SCF439/01

# Highlights

### Mimics baby's drinking rhythm




Babies know how to drink best! That's why our hands-free electric breast pump mimics their natural drinking rhythm - so you can find your own rhythm for effective\* expression and optimal milk output. The right size breast shields and shield inserts support your expression with a secure and comfortable fit every step of the way.

### Unique SkinSense shields



Our unique SkinSense silicone breast shields are designed for ultimate comfort, softly shaping around your breasts using your body's natural warmth. The breast shields are made with food-grade silicone for complete peace of mind.

### Transparent collection



The see-through collection cups provide an unobstructed view inside. They help you correctly place your nipples for your best milk

output, and visually track the volume of breast milk. For confident pumping from the first drop until finish.

### Ultra-light collection cup



One collection cup weighs just a little over 100 grams when empty - up to 3x lighter than a fully wearable breast pump. The cups sit comfortably in-bra with barely a sign they're there.

### Nipple sizing card



Use the nipple sizing card to find your ideal size of breast shields - and shield inserts if you use them. Measuring only takes a moment and with a range that fits 99% of moms, we have your comfortable and secure fit covered.

### Easy to clean and assemble




We designed our pump and accessories to be as simple as possible. With only a few connecting pieces, each collection cup is a breeze to take apart, clean (by hand or in a dishwasher), disinfect and put back together.

### Extra compatibility



These Hands-free Collection Cups can be used with the Philips Avent Premium Electric Breast Pump, Advanced Electric Breast Pump, and Hands-free Electric Breast Pump. Switching cups takes seconds, the benefits go on.

### Fits 99% of all moms



We support every mom to feel their best and having the right fit is an essential part of that. This starter set includes our full range of breast shields and shield insert sizes, fitting 99% of moms for ultimate comfort while pumping.

# Specifications

### Material

Cup, front cap, breast shield cover:

Polypropylene

Breast shield, membrane, valve, insert: Liquid silicone

Silicone tube: Silicone

### What is included

Collection cup (assembled): 2 pcs

Silicone tubes: 2 pcs

Breast shield (dust) cover: 2 pcs

Breast shield (21 mm, 24 mm): 2 pcs each

Insert 19 mm: 2 pcs

### Ease of use

**Breast pump use:** Easy to use, simple to clean and assemble

### Functions

**Flexibility in movement:** Hands-free pumping

**Breast shield and insert specifics:** Several sizes available for the best fit

