

# PHILIPS



IPL Hair removal  
device

Lumea IPL 7000  
Series

5 manual intensity settings

1 attachment: body

Lumea IPL App

Corded use



BRI920/00

## Great results, great value

Enjoy 12 months of hair-free smooth skin\*

Get hair-free smooth skin with our Lumea IPL 7000 Series. Treatments are gentle and effective with attachments for each body area.

### **Gentle and effective for long-lasting smooth skin**

- Fast results with treatments only every 2 weeks
- Developed with dermatologists to be easy and effective
- Extra-long cable for added flexibility while treating

### **Full solution for face and body**

- 1 attachment for body
- Comes with Satin Compact pen trimmer
- Suitable for a wide range of skin tones and hair colours

### **Personalised IPL hair removal treatment**

- Skin Tone Sensor and 5 intensity settings
- Optimise your routine with the Philips Lumea IPL app

## Lumea IPL 7000 Series

IPL Hair removal device

BRI920/00

# Highlights

### 2 x a month for fast results




Get started with the initial phase of 4 treatments only every 2 weeks — that's half as many treatments as other brands. Then touch up monthly to maintain the results.

### 1 attachment for body




The body attachment can be used for gentle and effective treatment. The wide treatment window is perfect for quickly treating larger body areas like the legs, arms and stomach.

### Philips Lumea IPL app




Our free coaching app helps you plan and stick to your treatment schedule, then takes you through each session step by step. Downloaded by more than 2.1 million users.

### Skin Tone Sensor



Choose from 5 intensity settings for a comfortable experience. The Skin Tone Sensor prevents you from treating areas of your skin that are too dark for IPL treatment.

### Developed with dermatologists




As a leader in health technology, Philips developed Lumea IPL in consultation with dermatologists, for easy and effective use in the safety of your home. Derived from technology used in professional salons, Lumea IPL offers gentle treatment, even on sensitive areas.

### Satin Compact pen trimmer




Complete your beauty routine: The Satin Compact pen trimmer is a discrete on-the-go beauty tool to quickly and easily remove even the finest facial hairs.

### Extra-long cable



Convenient to use thanks to an extra long cable for easy access and enhanced maneuverability.

### For many skin and hair types



IPL needs contrast between the pigment in the hair colour and the pigment in the skin tone; it therefore works on naturally dark blonde, brown and black hair and on skin tones from fair to medium brown (I-IV).

# Lumea IPL 7000 Series

IPL Hair removal device

BRI920/00

## Specifications

### Technical specs. attachments

**Body attachment:** Window size: 4.0 cm<sup>2</sup>, Treatment for body: legs, arms, stomach and underarms

### Safety and adjustable settings

**Integrated UV filter:** Protects skin from UV light

**5 light energy settings:** Adjustable to your individual skin tone

**Integrated safety system:** Prevents unintentional flashing

### Application mode

**Slide and Flash:** For quick application

**Stamp and Flash:** For treatment on small areas

**Corded/cordless use:** Corded use

### Pen trimmer

**Accessories:** Cleaning brush, Eyebrow comb

**Ease of use:** Compact handle

**Power:** AAA battery

### Power

**Corded use**

### Technical specifications

**High performance lamp:** Built to last, 250,000 flashes, equivalent to 20 years of lamp lifetime\*\*

**Voltage:** 100-240 V

### Attachments

**Body attachment (4 cm<sup>2</sup>):** For use below the neckline

### Service

**Warranty:** 2-year global warranty + 1 extra year of warranty upon product registration within 90 days

### Items included

**Instructions for use:** User manual

**Storage:** Beauty pouch

**Adapter:** 24 V/1500 mA

**Accessory:** Cleaning cloth

### Application time

**Bikini:** 4 min

**Lower leg:** 15 min

**Underarms:** 2.5 min

**Bikini line:** 4 min

**Face areas:** 2 min

