

PHILIPS

sonicare



Rechargeable
toothbrush

Series 6400

1000% more plaque removal*

Gentle on sensitive teeth and
gums

Next-Generation Sonicare
Technology

Visual Pressure Sensor



HX7411/04

Advanced care for sensitive teeth and gums

Gently removes 1000% more plaque*

Thorough cleaning can still be gentle on sensitive teeth and gums with 1000% more plaque removal. With our Next-Generation Sonicare Technology for a powerful yet gentle clean, even in those hard-to-reach areas.

Brushing made effortless

- Choose your ideal cleaning experience
- Brushing insights in the palm of your hand
- Keep brushing your best with brush head replacement reminder
- Guided brushing sessions

Designed for your lifestyle

- Carry the care, take it anywhere
- 21 days of regular brushing

Feel the care of Sonicare

- Extra gentle on sensitive teeth
- Next-Generation Sonicare Technology
- Help protect your gums with pressure sensor

Series 6400

Rechargeable toothbrush

HX7411/04

Highlights

Extra care for sensitive teeth



Thorough cleaning can still be gentle on sensitive teeth and gums. This S2 Sensitive brush head features a unique design with long, thin, ultra-soft bristles that keep every brush stroke kind. With more than 3,000 densely packed bristles, it helps you remove 1000% more plaque than a manual toothbrush.


Next-Generation Sonicare Technology



This electric toothbrush uses our Next-Generation Sonicare Technology which delivers our most reliable and consistent brushing for excellent results every time. The motor adjusts power automatically so there is zero drop in performance when you reach those harder-to-clean areas. Enjoy a gentle yet effective clean for your teeth and gums with 62,000 bristle

movements per minute. Sonicare Fluid Action supports the bristles to clean by driving fluid deep between teeth and along the gumline.

Visual pressure sensor




This Sonicare electric toothbrush has a light ring at its base that gently lets you know if you're overdoing the pressure. Simply ease off when it lights up and your gums will stay protected.

Modes & Intensities



Give your clean a boost with 9 brushing settings. Whether you want a little more vigor or a specific focus for your clean, you can have it. Choose between Clean, Sensitive and White mode. Pick one of three intensity settings.

Philips Sonicare app



Achieve your oral-health care goals with the Sonicare toothbrush and app. The two pair seamlessly to deliver guided brushing, tips and tricks, and content personalized to you. Together, you are unstoppable.

Travel case



Our travel case allows you to take your device on the go. It is strong enough to keep your Sonicare toothbrush safe, yet compact enough to slip into any bag, making it the perfect travel companion.

Series 6400

Rechargeable toothbrush

HX7411/04


Highlights

SmarTimer & BrushPacer



Two minutes is all you need for a complete clean, so our SmarTimer lets you know when time's up. For an even clean throughout your mouth, the BrushPacer divides your session into segments, with a slight pause in vibration being your sign to move on. This supports you in achieving an even clean every time you brush.

Replacement reminder



Did you know brush heads become less effective after three months? Dental professionals recommend changing your brush head regularly, that's why Philips Sonicare toothbrushes have brush head replacement reminder. It tracks how

often and how hard you brush, then reminds you to replace your brush head once it's time for a new one.

21-day battery life



Experience up to 21 days of regular brushing after a single charge. Bringing a new level of convenience to your brushing routine.

Series 6400

Rechargeable toothbrush

HX7411/04

Specifications

Intensities

High: To boost your clean
Medium: For an everyday clean
Low: For sensitive teeth and gums

Power

Voltage: 100-240 V

Technical specifications

Operating time (full to empty): 21 days
Battery: Rechargeable
Energy Consumption: Standby without display
<0.5 W
Battery type: Lithium ION

Design and finishing

Color: Black

Service

Warranty: 2-year limited warranty

Compatibility

Android compatibility: Android phones with Bluetooth 4.0 and above, Search Sonicare on Google Play Store for detailed OS compatibility.
Bluetooth wireless technology: Connected brushing app
iOS compatibility: iPhone with Bluetooth 4.0 and above, Search Sonicare on App Store for detailed iOS compatibility.

Ease of use

Handle: Slim ergonomic design
Handle compatibility: Easy click-on brush heads
Battery indicator: Illuminated icon indicates battery life
Timer: BrushPacer and SmartTimer
USB-A Charging: Power adapter not included

Items included

Handle: 1 6400 rechargeable toothbrush
Brush head: 1 S2 Sensitive
Travel case: 1 Travel case
Charger: 1 USB-A charger

Cleaning performance

Plaque removal: 1000% as effective*
Speed: 62,000 brush movements / min

Modes

Clean: For an exceptional everyday clean
Sensitive: For sensitive teeth and gums
White: To remove surface stains

Smart sensor technology

Pressure feedback: Light ring illuminates purple, Vibration and pulsing sound
Replacement Reminder: BrushSync lets you know when to replace your brush head

Sonicare app

Tracking and progress reports: Monitors brushing patterns over time, helping you maintain good oral hygiene.

