

## Lumify device compatibility lists

**Right. Where you need it.** Philips Lumify handheld ultrasound transforms your compatible Android or Apple smart device into a powerful, portable ultrasound system. With advanced processing built directly into the transducer, Lumify delivers exceptional image quality and rapid clinical insights—anytime, anywhere.

See compatible Android devices in the table below and compatible iOS devices in the table on the back.

### Lumify Android supported devices\*

Device name	Screen size	Connectivity options	Internal storage	Approximate continuous scan time	USB type	Supports Reacts	Supported OS
Getac ZX10	10 in	Wi-Fi / LTE	128 GB	4 h 39 min	C	Yes	12
Samsung Galaxy S24 Ultra	6.8 in	Wi-Fi / 5G	256 GB	4 h 51 min	C	Yes	14
Samsung Galaxy S24+	6.7 in	Wi-Fi / 5G	256 GB	6 h	C	Yes	14
Samsung Galaxy S24	6.1 in	Wi-Fi / 5G	256 GB	4 h 15 min	C	Yes	14
Samsung Galaxy S23	6.1 in	Wi-Fi / 5G	256 GB	4 h 41 min	C	Yes	14
Samsung Galaxy S9 Tab	11.0 in	Wi-Fi / 5G	128 GB	5 h 55 min	C	Yes	13
Samsung Galaxy S9+ Tab	12.4 in	Wi-Fi / 5G	256 GB	7 h 20 min	C	Yes	13
Samsung Galaxy S9 Ultra Tab	14.6 in	Wi-Fi / 5G	256 GB	7 h 40 min	C	Yes	13
Samsung Galaxy A9+	11.0 in	Wi-Fi / 5G	128 GB	5 h 5 min	C	Yes	13
Samsung Galaxy Tab Active5	8 in	Wi-Fi / 5G	128 GB	4 h 15 min	C	Yes	14
Samsung Galaxy Tab Active4 Pro (SM-T630)	10.1 in	Wi-Fi / 5G	64 GB	4 h 16 min	C	Yes	Out of the box
Samsung Galaxy Tab Active3	8 in	Wi-Fi / 5G	64 GB	4 h 14 min	C	Yes	Out of the box
Samsung Galaxy Tab A8	10.5 in	Wi-Fi / 5G	32 GB	3 h 57 min	C	Yes	Out of the box
Samsung A7 Lite	8.7 in	Wi-Fi	32 GB	3 h 20 min	C	Yes	Out of the box
Huawei Mate 30 Pro 5G	6.53 in	Wi-Fi / 5G	256 GB	3 h 57 min	C	Yes	Out of the box
Huawei Changxiang 2	10.1 in	Wi-Fi	128 GB	3 h 31 min	C	Yes	Out of the box
Panasonic Toughbook (FZ-S1)	7 in	Wi-Fi / LTE	64 GB	2 h 16 min	C	Yes	Out of the box
Panasonic Toughbook (FZ-A3)	10.1 in	Wi-Fi / LTE	128 GB	2 h 19 min	C	Yes	Out of the box
Samsung Galaxy Tab S2 8" (Only Model SM-T713)	8 in	Wi-Fi	32 GB	2 h 30 min	B	Yes	Out of the box
Samsung XCover Pro	6.3 in	Wi-Fi / LTE	64 GB	2 h 37 min	C	Yes	Out of the box

Items in **red** are new to this document revision.

Items in **blue** are compatible, but no longer available for purchase from the manufacturer.

\*Devices not included in the preceding table have not undergone compatibility testing with Lumify. Philips cannot guarantee Lumify system performance with non-compatible devices as compatibility testing with all devices on the market is not feasible. The devices have been evaluated using the operating system (OS) versions they were originally equipped with.

Subsequent OS releases, known as maintenance releases, that are higher than the version installed at the time of device launch have not been explicitly verified. Conducting explicit tests for newer major OS releases on older devices, which were initially equipped with a previous OS version, is not feasible. Lumify system functionality and performance may be impacted if significant operating system changes have been applied.

# Lumify iOS supported devices\*

Device name	Year released	Screen size	Connectivity options	Internal storage	Philips case available	LPM required	Connector type	Supported OS	Lumify SW version tested
iPad Air 11-inch M2	2024	11 in	Wi-Fi / LTE	128 GB	No	Yes	C	17.5.1	5.0
iPad Air 13-inch M2	2024	13 in	Wi-Fi / LTE	128 GB	No	Yes	C	17.5.1	5.0
iPad Pro 12.9-inch 6th Gen	2022	12.9 in	Wi-Fi / LTE	128 GB	No	Yes	C	17.3	5.0
iPad Pro 11-inch M4	2024	11 in	Wi-Fi / LTE	256 GB	No	Yes	C	17.5.1	5.0
iPad Air 5th Gen M1	2022	10.9 in	Wi-Fi / LTE	64 GB	No	Yes	C	17.2	5.0
iPad Pro 13-inch M4	2024	13 in	Wi-Fi / LTE	512 GB	No	Yes	C	17.5.1	5.0
iPhone 15	2023	6.12 in	Wi-Fi / LTE	128 GB	No	Yes	C	17.0.3	5.0
iPhone 15 Plus	2023	6.69 in	Wi-Fi / LTE	128 GB	No	Yes	C	17.0.3	5.0
iPhone 15 Pro	2023	6.12 in	Wi-Fi / LTE	128 GB	No	Yes	C	17.2.1	5.0
iPhone 14	2022	6.06 in	Wi-Fi / LTE	128 GB	No	Yes	Lightning	16.6	5.0
iPhone 14 Pro	2022	6.1 in	Wi-Fi / LTE	128 GB	No	Yes	Lightning	16.3.1	5.0
iPad 10th Gen	2022	10.9 in	Wi-Fi / LTE	64/256 GB	No	Yes	C	16.3.1	5.0
iPad Mini 6th Gen	2021	8.3 in	Wi-Fi / LTE	64/256 GB	No	Yes	C	16.6	5.0
iPhone 13 mini	2021	5.4 in	Wi-Fi / LTE	128/256/512 GB	No	Yes	Lightning	Out of the box	2.0
iPhone 13	2021	6.1 in	Wi-Fi / LTE	128/256/512 GB	No	Yes	Lightning	Out of the box	2.0
iPad 9th gen	2021	10.2 in	Wi-Fi / LTE	64/256 GB	Yes	Yes	Lightning	Out of the box	2.0
iPhone 12	2020	6.1 in	Wi-Fi / LTE	64/128/256 GB	No	Yes	Lightning	Out of the box	2.0

**Note: Lumify is not currently compatible with iOS version 18.x when used in conjunction with a USB-C port device.**

Content in **red** are new to this document revision.

\*Devices not included in the preceding table have not undergone compatibility testing with Lumify. Philips cannot guarantee Lumify system performance with non-compatible devices as compatibility testing with all devices on the market is not feasible. The devices have been evaluated using the operating system (OS) versions they were originally equipped with.

Subsequent OS releases, known as maintenance releases, that are higher than the version installed at the time of device launch have not been explicitly verified. Conducting explicit tests for newer major OS releases on older devices, which were initially equipped with a previous OS version, is not feasible. Lumify system functionality and performance may be impacted if significant operating system changes have been applied.

Please note that Reacts tele-ultrasound solution is not available on Lumify when using the above iOS devices. All iOS devices currently compatible require use of Lumify Power Module and utilize USB-C transducer cables and will continuously scan for 2 to 2.5 hrs.

Devices on the supported list satisfy the following criteria: 1. Passed Philips hardware testing, 2. Is readily available from the device manufacturer.

Philips cannot test every device available on the market, so this list represents a range of devices to satisfy a variety of user needs with respect to form factor, price, battery life and connectivity options.

