

PHILIPS

sonicare

Sonic electric
toothbrush

FlexCare Platinum

3 modes, 3 intensities

1 brush head

With pressure sensor



HX9111/21

Deep, gentle cleaning

10x more plaque removal and healthier gums*

With FlexCare Platinum, even the deepest clean is gentle enough to protect your teeth and gums. Thanks to the AdaptiveClean brush head and intuitive pressure sensor, you get maximum plaque removal and a gentle experience however you brush.

Proven to improve oral health

- Improves gum health in just 2 weeks
- Remove up to 10x more plaque with AdaptiveClean*
- Safe and gentle for your oral health needs

Designed around you

- Timers encourage thorough 2 minutes of brushing
- Long battery life and travel case
- Alerts you if you're pressing too hard

Provides a superior clean


- AdaptiveClean brush head flexes along your gum line
- Deep Clean mode tackles your trouble spots

Optimize your brushing

- Our unique technology brings you a powerful yet gentle clean
- Choose from 3 modes and 3 intensity settings


Highlights

Remove up to 10x more plaque




Click on the AdaptiveClean brush head to experience our deepest ever clean. You get 4x more surface contact** to remove up to 10x more plaque in hard-to-reach areas,* and a truly deep clean along the gum line and between teeth.

Look after those gums



With optimal cleaning from your FlexCare, gums will become healthier in 2 weeks. Removing up to 7x more plaque along the gum line than a manual toothbrush, you'll get your healthiest smile.

3 modes, 3 intensity settings



With FlexCare Platinum Connected, you experience a truly deep clean. Our 3 intensity settings let you boost your clean, while 3 modes

meet your specific needs: Clean mode – for superior daily cleaning, White – to remove surface stains, and Deep Clean for invigorating freshness.

Extra deep cleaning



Our Deep Clean mode means you can really focus on specific areas, paying them extra attention. We've combined our unique brush head motion with an extended 3 minutes of brushing time so your teeth get a truly deep clean.

Too much pressure?




You might not notice if you're brushing too hard, but your FlexCare Platinum will. If you use too much pressure, your handle will gently pulsate. This is a reminder to ease off the pressure, and let your brush head do the work. 7 out of 10 people found this feature helped them become a better brusher.

Philips Sonicare technology



Powerful sonic vibrations whip up your toothpaste into plaque fighting bubbles, and drive them deep between your teeth and along your gum line. At the same time your teeth will experience 62,000 gentle but effective brush movements. You'll get a month's worth of manual brushing in just 2 minutes.

Up to 2x healthier teeth



Click on our AdaptiveClean brush head for a deep, gentle clean that leaves gums up to 2x healthier than a manual toothbrush. Soft rubber sides flex along every contour of your teeth and gums to absorb any excess pressure. Bristles can then reach deep between teeth giving you a comfortable, effective deep clean.


FlexCare Platinum

Sonic electric toothbrush

HX9111/21

Highlights

QuadPacer and Smartimer



It takes just 2 minutes to thoroughly brush your teeth. Our QuadPacer lets you know when you've spent the optimal amount of time on each section of your mouth while our Smartimer signals when your total time is up. Together they help you reach your recommended brushing time, every time.

Brush with confidence



Whatever your specific needs, you can be sure of a safe brushing experience with FlexCare Platinum. Our sonic technology is suitable for use with braces, fillings, crowns and veneers, and can even be used in the daily treatment of gum disease.

Designed to fit your lifestyle



Our premium travel case lets you store your toothbrush hygienically, while our compact travel charger keeps you topped up when you're on the go. You'll enjoy two weeks of regular use from a single full charge, so being on the move means you needn't miss out on that fresh feeling.

Specifications

Power

Voltage: 110-220 V

Technical specifications

Battery: Rechargeable

Operating time (full to empty): Up to 3 weeks**

Battery type: Lithium ION

Design and finishing

Color: Slate gray

Service

Warranty: 2-year limited warranty

Ease of use

Display: Illuminated display

Handle: Slim ergonomic design

Battery indicator: Light shows battery status

Brush head system: Easy click-on brush heads

Brushing time: Up to 3 weeks***

Items included

Handle: 1 FlexCare Platinum

Brush heads: 1 AdaptiveClean

Travel case: 1

Charger: 1

Cleaning performance

Health benefits: Improves gum health in only two weeks

Performance: Removes up to 10x more plaque*

Pressure sensor: Vibrates handle to alert user

Speed: Up to 62000 brush movement/min

Whitening benefits: Helps naturally remove stains

Timer: Quadpacer and SmarTimer

Modes

Clean: For exceptional everyday clean

3 intensities: Low, Medium, High

White: Removes surface stains

Deep Clean: For an invigorating deep clean

