

PHILIPS

AVENT



Essential

Single Electric Breast Pump



SCD323/20

Express effortlessly

The Philips Avent electric breast pump keeps milk flow at an optimum and gently adapts to your size and shape. Quick, comfortable, and convenient for each new mother.

Personalized

- Personalised experience, 4 + 4 setting levels
- Minimal parts and intuitive setup
- Memory function
- Small and Compact
- Healthy storage of breast milk

Gentle and comfortable

- Soft and adaptive silicone cushion


Highlights

Personalised experience



Finely tune every session to your needs with a wide range of stimulation and expression settings. Our breast pump offers 4 stimulation and 4 expression levels for a personalised experience.

Memory setting




Know what you like? The breast pump automatically remembers your last settings, so all you need to do is sit down and press start.

Soft adaptive silicone cushion




One size fits all. The silicone cushion gently flexes and adapts to fit your nipple.

Small and Compact




Small and Compact motor unit

Easy to clean and set up



Our closed expression system means milk stays out of the tubing so there's less to clean. Fewer parts also make it a breeze to put back together again.

Healthy storage of breast milk



Convenient, pre-sterilized, one-time use bags provide safe and secure storage for a nursing mother's precious breast milk.

Specifications

Power

USB-C connection

Main voltage: 5V and ≥ 1 A

Supply unit is not included. Not chargeable, corded use only.

Material

0% BPA*

What is included

Milk bags: 5 pcs

Motor unit

Silicone tube

USB-C cable

Natural bottle (125 ml)

Expression kit

Screw-on lid

Closing plate: 1 pcs

Slow teat: 1 pcs

Breastpads: 4 pcs

Ease of use

Breast pump use: Easy cleaning & assembly, memory function, portable design

Functions

Settings: 4 expression levels, 4 massage levels

