

PHILIPS

Front cap

Hands-free

Breast pumps

Pink

CP2216



To connect the collection cups and motor unit

Check for compatibility below

This front cap secures the connection between the Philips Avent Hands-free electric breast pump collection cups and the motor unit. The cap can also be rotated between four different positions for comfortable pumping and breastfeeding.

Specifications

Replaceable part

Fits product types: SCD553, SCF439, SCF442,
SCF481, SCF531, SCF532, SCF533, SCF534,
SCF537, SCF539, SCF543, SCF547

