

# PHILIPS

## Silicone tube

Hands-free

Breast pumps

Transparent

CP2211



## Connecting the collection cups with the motor unit

Check for compatibility below

The silicone tube connects the Philips Avent Hands-free electric collection cups with the motor unit. Do not trim the tubes as this increases the suction level, which may hurt or feel uncomfortable and lower your milk production.

# Specifications

## Replaceable part

Fits product types: SCD553, SCF439, SCF442,  
SCF481, SCF531, SCF532, SCF533, SCF534,  
SCF537, SCF539, SCF543, SCF547

