

PHILIPS

sonicare



Rechargeable Toothbrush

DiamondClean Prestige
9900

SenseIQ

20 x more plaque removal*

15 x healthier gums**

Visual pressure sensor



HX9992/43

Personalised care for your oral health

Gently removes up to 20 x more plaque*

Experience our most advanced electric toothbrush: Philips Sonicare DiamondClean Prestige 9900 with SenseIQ technology. It enhances our proven sonic technology with intelligence, sensing and adapting to deliver personalised oral care.

Brushing made effortless

- Help protect your gums with our pressure sensor
- Choose your ideal cleaning experience
- Personalised guidance and real-time feedback

Designed for your lifestyle

- Designed to care, with care
- Keep brushing your best with the brush head replacement reminder
- 14 day battery life

Feel the care of Sonicare

- 20 x more effective on plaque*, gentle on gums
- SenseIQ Technology for your personalised experience
- Sonicare Fluid Action

DiamondClean Prestige 9900

Rechargeable Toothbrush

HX9992/43


Highlights

Our best plaque removal



Everyone cleans differently, so we designed this brush head with multi-angle bristles to keep plaque removal on course no matter your technique. Our All-in-One A3 brush head delivers our best cleaning and removes up to 20 x more plaque, 15 x healthier gums in six weeks and up to 100% more stain removal in less than two days versus a manual toothbrush.

SenseIQ Technology




SenseIQ Technology helps balance your brushing style and senses your pressure, motion and coverage up to 100 times per second while you brush. If you apply too much pressure, the toothbrush will step in and automatically adjust the intensity level, helping to protect your teeth and gums.

Sonicare Fluid Action




Philips Sonicare toothbrushes clean gently yet effectively and care for your teeth and gums with 62,000 bristle movements per minute. Sonicare Fluid Action supports the bristles to clean by driving fluid deep between teeth and along the gumline.

Visual pressure sensor



This Sonicare electric toothbrush has a light ring at its base that gently lets you know if you're overdoing the pressure. Simply ease off when it lights up to help your gums stay protected.

Modes and intensities



Enhance your oral care routine with 15 brushing settings. Whether you're looking for a deep clean or a specific focus, you can have it with this electric toothbrush. Select from five modes—Clean, White+, Gum Health, Deep Clean and Sensitive—and adjust to one of three intensity levels for a personalised clean. Simply open the app, select your preferences and let your toothbrush do the rest.

Sonicare App



Achieve your oral health goals with your Sonicare toothbrush and app. The two pair seamlessly to provide real-time coaching, personalised guidance and tailored insights to improve your brushing routine. Track your progress, refine your technique and receive continuous feedback to make the most out of each brushing session. Together, you are


unstoppable. Meanwhile, auto-syncing keeps brushing data up to date in the app, even when you don't brush with the app by your side.

Premium design and accessories




From its sleek design to its state-of-the-art charging stand, every detail of this electric toothbrush is crafted with attention to detail by our dedicated designers and engineers. Plus, with the included charging travel case, you can maintain your oral care routine, even while on the go.

Replacement Reminder



Did you know brush heads become less effective after three months? Dental professionals recommend changing your brush head regularly, so that's why Philips Sonicare toothbrushes have a brush head replacement reminder. It tracks how often and how hard you brush, then reminds you to replace your brush head once it's time for a new one.

Battery life



Experience up to 14 days of regular brushing on a full charge, bringing a new level of convenience to your daily routine.

DiamondClean Prestige 9900

Rechargeable Toothbrush

HX9992/43

Specifications

Intensities

High: To boost your clean

Medium: For an everyday clean

Low: For sensitive teeth and gums

Power

Power supply: 100-240 V

Technical specifications

Battery: Rechargeable

Operating time (full to empty): 14 days

Energy consumption: Standby without display
<0.06 W

Battery type: Lithium ION

Design and finishing

Colour: Black

Service

Warranty: 2-year limited warranty

Software support: Philips offers relevant software updates for a period of 2 years after the date of purchase.

Compatibility

Android compatibility: Android phones with Bluetooth 4.0 and above

Bluetooth wireless technology: Connected
Sonicare app

iOS compatibility: iPhone with Bluetooth 4.0 and above

Ease of use

Handle: Sleek and compact design

Battery indicator: Illuminated icon indicates battery life

Handle compatibility: Easy click-on brush heads

Timer: BrushPacer and Smartimer

Items included

Handle: 1 DiamondClean Prestige 9900

Brush heads: 4 A3 All-in-One

Travel case: 1 Prestige USB charging travel case

Charger: 1 Charging base and stand

Cleaning performance

Plaque removal: 20 x more effective*

Gum Health: Up to 15 x healthier gums***

Whitening: Up to 100% more stain removal***

Speed: 62,000 brush movements/min

Modes

Clean: For an exceptional everyday clean

Deep Clean: For an invigorating deep clean

Gum Health: Gently massages your gums

Sensitive: For sensitive teeth and gums

White+: To help remove surface stains

Smart sensor technology

Pressure feedback: Ring lights up purple, Vibration and pulsing sound

SenseIQ: Automatically adjusts brushing based on pressure and motion

Replacement Reminder: Lets you know when to replace your brush head

Sonicare app

Tracking and progress reports: Track your brushing habits and improvements over time with in-depth data and personalised insights.

Real-time coaching with SenseIQ: Receive instant feedback on pressure, motion and coverage to improve your brushing technique as you go.

Personalised guidance: Get customised brushing tips based on your brushing patterns and habits.

