

PHILIPS

sonicare



Sonic electric
toothbrush

4100 Series

Sonicare technology
QuadPacer and SmarTimer
Slim ergonomic design
14 day battery life



HX3689/42

Goodbye manual toothbrush. Hello Sonicare.

Helps whiten teeth every time you brush

Sonic technology combined with our brush action gently removes stains for naturally whiter teeth every time you brush.

Designed around you

- Use less. Reuse more
- Switch between 2 intensity settings
- Always know when to replace your brush heads
- 14 day long-lasting battery life

A whiter, healthier smile

- Our unique technology brings you a powerful yet gentle clean
- Helps whiten teeth every time you brush

Optimise your brushing

- Optimise your brushing with SmarTimer and QuadPacer

Protects your gums

- Pressure Sensor helps protect your teeth and gums


4100 Series

Sonic electric toothbrush

HX3689/42


Highlights

Helps whiten teeth every time you brush



Sonicare electric toothbrush with advanced sonic technology gently removes stains for naturally whiter teeth every time you brush.

Protects your gums



Built-in pressure sensor automatically detects the pressure you apply, warns you and reduces the toothbrush vibrations automatically to help protect your gums. The toothbrush will make a pulsing sound as a reminder to ease off the pressure. 7 out of 10 people found this feature helped them become better brushers.

SmartTimer and QuadPacer




2 minute SmartTimer and 30 second QuadPacer guide you to brush for the recommended time in all areas of your mouth to ensure a complete clean.

2 intensity settings



Switch between 2 intensity settings to make the brushing session comfortable for you.

Advanced Sonicare Technology



Powerful bristle vibrations drive microbubbles deep between your teeth and along your gum line for a refreshing experience. You will get two months' worth of manual brushing in just 2 minutes.* 31,000 brush strokes per minute gently clean your teeth, break up plaque and sweep it away for an exceptional daily clean.

BrushSync replacement reminder



All brush heads wear out over time, so you'll want to keep an eye on yours to make sure you're still getting a great clean. Our BrushSync technology tracks how long you've been using your brush head for and how hard you've been

brushing. The BrushSync replacement reminder on your handle and a short beep remind you when it's time to replace it.

Long-lasting battery life



Battery life lasts up to 14 days so you can go long periods without charging.

Pulling the plug on waste



Our ambition at Philips is to reduce electronic waste. One way we do this is to minimise the number of adapters we put in the market. These products come without a power adapter because we are committed to reducing electronic waste.

4100 Series

Sonic electric toothbrush

HX3689/42

Specifications

Power

Voltage: DC5V

Technical specifications

Battery: Rechargeable

Operating time (full to empty): 14 days**

Energy consumption: Standby without display
<0.5 W

Battery type: Lithium ION

Design and finishing

Colour: Malibu Prismatic

Service

Warranty: 2-year limited warranty

Ease of use

Handle: Slim ergonomic design

Brush head system: Easy click-on brush heads

Battery indicator: Light shows battery status

BrushSync Replacement Reminder: Replacement reminder, Always know when to, replace brush heads, Reminder icon lights up, to always ensure best results

Items included

Handle: 1 4100 Series

Brush heads: 1 W Optimal White

Charger: 1

Cleaning performance

Performance: Helps whiten teeth every time you brush

Timer: Quadpacer and Smartimer

Pressure feedback: Vibrates handle to alert user

Modes

Clean: For exceptional everyday cleaning

2 intensities: High, Low

Smart sensor technology

BrushSync Replacement Reminder: Always know when to, replace brush heads

Pressure sensor: Alerts when brushing too hard

