

PHILIPS

sonicare



Rechargeable
toothbrush

1100

Gently removes 3 x more plaque

Easy Start

SmartTimer and QuadPacer

14 day battery life



HX3901/03

Essential Care for healthy teeth and gums

Gently removes up to 3 x more plaque*

Switch to electric brushing effortlessly with the essential Philips Sonicare 1000 Series. Enjoy up to 3 x more plaque removal vs a manual toothbrush

Brushing made effortless

- Easy switch to electric brushing
- Guided brushing sessions

Designed for your lifestyle

- 14 day battery life

Feel the care of Sonicare

- Gently removes up to 3 x more plaque*
- Sonicare Fluid Action

1100

Rechargeable toothbrush

HX3901/03

Highlights

Say goodbye to plaque




This electric toothbrush comes with our InterCare brush head. The extra-long bristles help remove more plaque from between teeth and in hard-to-reach spots.

Sonicare Fluid Action



Philips Sonicare toothbrushes clean gently yet effectively and care for your teeth and gums with 31,000 bristle movements. Sonicare Fluid Action supports the bristles cleaning by driving fluid deep between teeth and along the gumline.

Easy Start



Sonicare easy-start programme is designed for those who are new to brushing with an electric toothbrush. The programme gives you the option of a gradual, gentle increase in brushing power across your first 14 cleans.

SmartTimer and QuadPacer



Ace your session time with Sonicare's cleaning timers. Every 30 seconds, the BrushPacer will prompt you to clean a new area. After 2 minutes, the SmartTimer will indicate your session is complete.

Battery life



Experience up to 14 days of regular brushing after one full charge. Bringing a new level of convenience to your brushing routine.

1100

Rechargeable toothbrush

HX3901/03

Specifications

Intensities

Medium: For an everyday clean

Power

Voltage: 100-240 V

Technical specifications

Battery: Rechargeable

Operating time (full to empty): 14 days

Battery type: Lithium ION

Energy consumption: Standby without display
<0.5 W (automatically reached within 1 minute)

Design and finishing

Colour: Marine Blue

Service

Warranty: 2-year limited warranty

Ease of use

Handle: Slim ergonomic design

Handle compatibility: Easy click-on brush heads

Battery indicator: Illuminated icon indicates
battery life

Timer: SmarTimer and QuadPacer

Items included

Handle: 1 1100 rechargeable toothbrush

Brush head: 1 InterCare

Charger: 1 Charging cable

Cleaning performance

Plaque removal: Up to 3 x more plaque removal*

Modes

Clean: For exceptional everyday cleaning

Consumer Trade Item

Height: 22 cm

Width: 8 cm

Length: 4.5 cm

Net weight: 0.369 kg

Gross weight: 0.421 kg

GTIN: 8720689043762

Country of origin: CN

Outer Carton

Height: 25 cm

Width: 10 cm

Length: 15.8 cm

Net weight: 1.107 kg

Gross weight: 1.375 kg

GTIN: 18720689043769

Pallet Layer

Height: 25 cm

Width: 80 cm

Length: 120 cm

Pallet

Height: 112.8 cm

Width: 80 cm

Length: 120 cm

Net weight: 265.68 kg

Gross weight: 341.26 kg

GTIN: 68720689043764

