

# PHILIPS

## sonicare



2-pack brush heads

I InterCare

Reaches deep in between teeth

Medium soft



HX9002/87

## The brush head for everyday cleaning

### Reaches deep in between teeth

Target hidden plaque caught between teeth with this InterCare replacement brush head from Philips Sonicare. For a great clean every day.

#### Care for the planet

- Made with 70% bio-based plastic
- Paper-based recyclable packaging

#### Easy-to-use in every way

- Fits all Philips Sonicare click-on handles
- Keep brushing your best with BrushSync reminders

#### For optimal performance

- Deep cleaning between teeth and in hard-to-reach areas
- Proven to deliver effective oral care
- Sonicare Fluid Action

# Highlights

## Clean between teeth




Our InterCare brush head uses extra-long bristles to help remove up to 3 x more plaque\*\* from in between teeth and in hard-to-reach spots.

## Sonicare Fluid Action




Philips Sonicare toothbrushes gently yet effectively clean and care for your teeth and gums with up to 62,000 bristle movements. Sonicare Fluid Action supports the bristles to clean by driving fluid deep between teeth and along the gumline.

## Clinically tested



All Philips Sonicare brush heads are meticulously designed and clinically tested. So you can be sure that they meet the highest standard in quality, safety and performance.

## Replacement reminders



Did you know brush heads become less effective after three months? Dental professionals recommend changing your brush head regularly, that's why Philips Sonicare toothbrushes have BrushSync technology. It tracks how often and how hard you brush, then reminds you to replace your brush head once it's time for a new one.

## Click-on design




This brush head is compatible with all Philips Sonicare electric toothbrushes.\*\*\* Simply click it on or off for easy replacement and cleaning.

## Bio-based plastic



We are working to reduce our use of fossil-fuel-based plastic in our products. That's why 70% of the plastic in this brush head is bio-based.\*

## Paper-based packaging



All of our brush heads come in paper-based packaging which can be recycled where facilities are available.

# Specifications

## Design and finishing

**Colour:** White

**Bristle stiffness feel:** Medium soft

**Reminder bristles:** Bristles' blue colour fades away

**Size:** Standard

**Smart Brush Head Recognition**

## Compatibility

**Brush head system:** Click-on

**Not suitable for:** Philips One

**Compatible Toothbrushes:** Philips Sonicare

## Items included

**Brush heads:** 2 I InterCare

## Quality and performance

**Replacement:** Every 3 months

**Tested:** for optimal usage

**Benefit:** Plaque removal

