

PHILIPS

sonicare



Oral Irrigator

Compact Flosser 1000

Effective flossing in just 60 seconds*



HX3333/23

Effective flossing in just 60 seconds*

Removes up to 99.9% of plaque in treated areas*

The Philips Sonicare Compact Flosser is fast, effective, and ready to take anywhere. The flosser removes up to 99.9% of plaque* from between your teeth and along your gumline in just 60 seconds.

Compact for on-the-go

- Collapsible 200ml reservoir for easy storage and travel
- 14 days of battery life between charges
- 2-in-1 nozzle-storage and filling port

The easy way to start flossing

- Enjoy Pulse, Clean or Sensitive mode
- Two nozzles for gentle and effective cleaning
- Clean every quadrant, never miss a spot
- IPX7 water-resistant
- Two ways you can add water

Compact Flosser 1000

Oral Irrigator

HX3333/23

Highlights

Compact and portable



The flosser's compact design has a collapsible reservoir, making it easy to put away or pack on your next trip. The reservoir holds upto 200ml of water, enough for one complete 60-second clean.

Three flossing modes



Tailor any session to your needs with three flossing modes. There's Pulse to massage and stimulate gums. Clean for everyday use. And Sensitive for an extra touch of gentleness in your routine.


Standard and Comfort nozzles



The flosser comes with two nozzles to remove up to 99.9% of plaque in treated areas*. The N1 Standard nozzle gives a powerful clean between

teeth. While the N2 Comfort nozzle, with its soft tip, has been designed for sensitive gums. Each nozzle can rotate 360 degrees to help you floss hard-to-reach areas, like at the back of your mouth.

15-second pacer



With its 15-second pacer, the flosser will briefly pause waterflow to let you know when its time to move to the next quadrant. It gives your whole mouth the right amount of flossing time, every time.

Rechargeable



With a charge that lasts up to 14 days, the Compact Flosser is super handy on-the-go. You can easily charge it at home before a trip, or top it up while you're away with the USB-A charging cable.

Water-resistant



The Compact Flosser has an IPX7 water-resistance rating, making it truly bathroom friendly.

Nozzle storage



Make the most of the clever design by tucking a nozzle into the channel at the top of the flosser when not in use. The channel doubles as a filling port.

Two filling methods



You can fill the water reservoir easily in two ways. Either remove the reservoir and fill it up like a glass - ideal when it's empty. Or top it up using the filling port near the nozzle attachment.

Compact Flosser 1000

Oral Irrigator

HX3333/23

Specifications

Technical specifications

Energy Consumption: Off mode <0.3 W, Standby without display <0.5 W

Design and finishing

Color: Purple

Service

Warranty: 2-year limited warranty

Ease of use

Handle: Collapsible design

Nozzle attachment: Easily clicks on and off

Nozzle storage: Easy to store nozzle on-the-go

Performance

Cleaning: Clean entire mouth in 60 sec

Items included

USB charging cable: 1

Compact Flosser: 1

N1 Standard nozzle: 1

N2 Comfort nozzle: 1

Modes

Clean: For everyday clean

Sensitive: For a gentle clean

Pulse: For gum stimulation

Health benefits

Plaque removal: Removes up to 99% of plaque in treated areas

Consumer Trade Item

Height: 17.8 cm

Width: 13.2 cm

Length: 8.75 cm

Net weight: 0.326 kg

Gross weight: 0.485 kg

GTIN: 8720689046985

Country of origin: VN

Harmonized system code: 33069000

Outer Carton

Height: 15.5 cm

Width: 19.1 cm

Length: 19.3 cm

Net weight: 0.652 kg

Gross weight: 1.179 kg

GTIN: 18720689046982

Pallet Layer

Height: 15.5 cm

Width: 121.9 cm

Length: 101.6 cm

Pallet

Height: 105 cm

Width: 121.9 cm

Length: 101.6 cm

Net weight: 117.36 kg

Gross weight: 242.48 kg

GTIN: 68720689046987

