

PHILIPS

sonicare



Rechargeable
toothbrush - 2 pack

Advanced Clean

1000% more plaque removal*

Up to 500% healthier gums**

Next-Generation Sonicare
Technology

Pressure sensor



HX7129/01

Advanced care for healthier gums**

1000% more plaque removal*

Support your gum health every time you brush. Our Next-Generation Sonicare Technology offers a powerful yet gentle clean, even in hard-to-reach areas, with up to 500% healthier gums** and 1000% more plaque removal*.

Brushing made effortless

- Choose your ideal cleaning experience
- Guided brushing sessions
- Keep brushing your best with brush head replacement reminder

Designed for your lifestyle

- Carry the care, take it anywhere
- 21 days of regular brushing

Feel the care of Sonicare

- Up to 500% healthier gums**
- Next-Generation Sonicare Technology
- Extra support to protect your gums

Advanced Clean

Rechargeable toothbrush - 2 pack

HX7129/01

Highlights

Extra care for your gums



Support your gum health with the GC Gentle Clean brush head. The bristle layout allows the outer bristles to keep your gum line brushed even while you're focusing on your teeth. The result is up to 500% healthier gums in six weeks**. Its engineered flexibility also helps you remove up to 1000% more plaque than a manual toothbrush*.

Next-Gen Sonicare Technology



This electric toothbrush uses our Next-Generation Sonicare Technology which delivers our most reliable and consistent brushing for excellent results every time. Brushing power automatically adapts to deliver a gentle yet effective clean, even in hard-to-reach areas. Designed for consistency and durability Enjoy a gentle yet effective clean for your teeth and gums with 62,000 bristle movements per minute. Sonicare Fluid Action supports the bristles to clean by driving fluid deep between teeth and along the gumline.

Pressure sensor



It's easy to brush too hard, so this Philips Sonicare toothbrush has a smart optic sensor to detect excessive pressure. The sensor lets you know if you're overdoing it through haptic vibrations. Ease off and your gums will stay protected!

Modes & Intensities



Give your clean a boost with 4 brushing settings. Whether you want a little more vigor or a specific focus for your clean, you can have it. Choose between Clean and Gum mode. Pick one of two intensity settings.

Travel case




Our travel case allows you to take your device on the go. It is strong enough to keep your Sonicare toothbrush safe, yet compact enough to slip into any bag, making it the perfect travel companion.

SmartTimer & BrushPacer



Ace your session time with Sonicare's cleaning timers. Every 20 seconds, the BrushPacer will prompt you to clean a new area. After 2 minutes, the SmartTimer will indicate your session is complete.

Replacement reminder



Did you know brush heads become less effective after three months? Dental professionals recommend changing your brush head regularly, that's why Philips Sonicare toothbrushes have brush head replacement reminder. It tracks how often and how hard you brush, then reminds you to replace your brush head once it's time for a new one.

21-day battery life



Experience up to 21 days of regular brushing after a single charge. Bringing a new level of convenience to your brushing routine.

Advanced Clean

Rechargeable toothbrush - 2 pack

HX7129/01

Specifications

Intensities

High: To boost your clean

Low: For sensitive teeth and gums

Power

Voltage: 100-240 V

Technical specifications

Operating time (full to empty): 21 days

Battery: Rechargeable

Energy Consumption: Standby without display
0.08 W

Battery type: Lithium ION

Design and finishing

Color: White

Service

Warranty: 2-year limited warranty

Ease of use

Handle: Slim ergonomic design

Handle compatibility: Easy click-on brush heads

Battery indicator: Illuminated icon indicates
battery life

Travel case: Travel case

Timer: BrushPacer and SmartTimer

USB-A Charging: Power adapter not included

Items included

Handle: 2 Advanced Clean rechargeable
toothbrush

Brush heads: 3 GC Gentle Clean

Travel case: 2 Travel case

Charger: 2 USB-A chargers

Cleaning performance

Plaque removal: 1000% more effective*

Gum Health: Up to 500% healthier gums**

Speed: 62,000 brush movements / min

Modes

Clean: For exceptional everyday clean

Gum Health: Gently massages your gums

Smart sensor technology

Pressure feedback: Light ring illuminates purple,
Vibration and pulsing sound

Replacement Reminder: BrushSync lets you know
when to replace your brush head

