

PHILIPS

AVENT

Bottle Warmer Gift Set

Premium

Warms evenly, no hotspots

Quick warming and defrosting mode

Suitable for milk and baby food

SCF358/10



Smart and easy warming

Take the guesswork out of feed preparation

Prepare warmed feeds in just 3 minutes* with a bottle warmer that manages the temperature for you. The smart temperature control sensor automatically adjusts the heating pattern so the warmer warms quickly and evenly.

Easy to clean

- Easy to clean

Automatically warms to the ideal temperature

- Keeps milk warm for up to 60 min and automatically shuts off

Defrosts milk and baby food containers

- Warms up baby food containers as well as milk

AirFree vent designed for less air ingestion

- Natural Response Bottle with AirFree Vent

For Philips Avent bottles and other leading brands

- Compatible with most leading bottle and baby food jar brands

Highlights

Manages temperature for you




Set the milk volume, press start and let the smart temperature control take care of the rest for you. It detects the milk's initial temperature and quickly warms it to the ideal temperature and keeps it there for up to 60 minutes.

Uses hospital-standard water-bath warming



Warm milk the way nurses do in hospitals using water-bath warming. Compared to other methods, it better preserves milk proteins essential for building up a baby's immune system. Gentle and quick, our technology helps promote your baby's healthy development with every feed.

Ready whenever you are



Our bottle warmer keeps milk warm for up to 60 minutes, just in case you need more flexibility around mealtime. Then automatically shuts off, giving you peace of mind.

Easily defrosts frozen feeds



Like to keep the freezer stocked with extra feeds? The bottle warmer quickly defrosts milk and baby food containers too.

Supports transition to solids



When your little one is ready to graduate to solids, the bottle warmer also defrosts and warms baby food containers.

Less time in the kitchen



Designed with just one piece so cleaning is easy and you can enjoy more time with your little one.

Fits your baby's favourites



Designed to fit your baby's favourite Philips Avent bottles and most leading baby bottle and baby food jar brands.

Natural Response Bottle with AirFree Vent



The Natural Response Teat works with your baby's natural feeding rhythm**: the teat has a unique opening which only releases milk when baby actively drinks. The AirFree vent is designed for extra protection against colic, gas and reflux by preventing air from getting into your little ones tummy while upright feeding.

Specifications

Product Material

Bottle warmer: PP, ABS

Baby bottle: PP

Teat: Silicone

Technical specifications

Power consumption (standby mode): <0.3 W
(period reaches automatically standby mode:
<1 min)

Power consumption: 400 W

Voltage: 220-240 V ~ 50-60 Hz

Weight and dimensions

Product dimensions (W x H x D):

133 x 160 x 148 mm

Retail pack dimensions (WxHxD):

240 x 185 x 148 mm

Country of origin

Bottle warmer: Made in China

Baby bottle: Made in Indonesia

What is included

9 oz/260 ml Baby bottle: 1 pcs

AirFree vent: 1 pcs

Medium Teat: 1 pcs

Development stages

Baby bottle: 3-6 m

Packaging specifications

Paper-based packaging*

