

# PHILIPS

## Knee Massager

Vibration Massage

Heated Knee Wrap

Touch Screen Control

Adjustable Velcro

PPM5303



## Vibrating Heat , Knee Relaxation

Take knee issues seriously

High-frequency vibration , heating relaxation Warm soothing relief, warms knees

### **High-Frequency Vibration**

- 3 vibration modes
- High-Frequency Vibration Relieves Knee Pressure
- Warm Soothing Relief Warms Knees

# Highlights

## 3 vibration modes

Auto, Low, High 3 levels massage intensity

## High-Frequency Vibration Relieves Knee Pressure

Daily Massage Soothing Relaxation Deep Massage

## Warm Soothing Relief Warms Knees

3 Adjustable Heat Levels Low: 45°C / Medium: 50°C / High: 52°C\*

# Specifications

## General specifications

Net weight: 0.97kg

Color: Blue

Product dimension: 243x162x145mm mm

Body part to use: Both knees

Functional Components: 2 vibration motor, heating film

## Main parameters

Temperature of heating: 50/52 °C

Battery: 3000mAh

Charging port: type-C

Charging Voltage: 5V

Charging Current: 2A

Operating Power: 10W

Operating Voltage: DC 3.7V

## Country of origin

Made in: China

## Service

Warranty period: 2 years

