

# PHILIPS

## Neck massager

### Neck & Shoulder massager

Bidirectional Massage

Heat for Deeper Relief

Simple Operation

Durable Material

PPM3528BK



## Neck & Shoulder Massager with Hot Compress

The PPM3528 is a neck and shoulder massager with heat function, equipped with speed-adjustable and automatically rotating massage heads, delivering a realistic massage experience that mimics human hands.

### Suitable for multiple body areas

- Can be used on the shoulder, neck, back and legs

### Soothe aching muscles

- Soothing Heat for Deeper Relief

### Simple Operation

- Power Button, Mode Button, Intensity Button, Heating Button

### Drive away fatigue

- Deeply kneading to drive away fatigue

### Durable Material

- Made with breathable mesh fabric and soft PU leather

# Neck & Shoulder massager

Neck massager

PPM3528BK/27

## Highlights

### Four Buttons for Simple Operation

1.Power Button: Press and hold for 1.5 seconds to turn on/off 2.Mode Button: Switch between 3 massage modes 3.Intensity Button: 3 adjustable speed levels 4.Heating Button: 2 adjustable heating temperatures\* 40°C±3°C / 45°C±3°C

### Bidirectional Massage

Equipped with adjustable speed, automatic rotation, provide life-like massage experience.

### Skin-friendly & Durable Material

Made with breathable mesh fabric and soft PU leather

### Multifunctional Massager

Can be used on the shoulder, neck, back and legs

### Heat for Deeper Relief

Soothe aching muscles and simultaneously improve blood circulation

