

PHILIPS

Knee massager

Vibration Massage

Heated Knee Wrap

Touch Screen Control

Adjustable Velcro

PPM5303



Vibrating Heat Knee Relaxation

Take knee issues seriously

High-Frequency Vibration Relieves Knee Pressure Warm Soothing Relief Warms Knees

High-Frequency Vibration Relieves Knee Pressure

- Warm Soothing Relief Warms Knees
- High-Frequency Vibration Relieves Knee Pressure
- Touch Screen Easy Operation

Highlights

Warm Soothing Relief Warms Knees

3 Adjustable Heat Levels Low: 45°C/Medium: 50°C/High: 52°C*

High-Frequency Vibration Relieves Knee Pressure

3 level intensity adjustment, Auto → Low → High

Touch Screen Easy Operation

3 Sleek Buttons Temperature Indicator Intensity Indicator

Specifications

General specifications

Net weight: 0.97 kg

Product dimensions: 243 x 162 x 145 mm mm

Colour: Grey

Body part to use: Both knees

Functional Components: 2 vibration motors, heating film

Main parameters

Temperature of heating: 50/52 °C

Battery: 3000 mAh

Charging port: type-C

Charging Voltage: DC 5 V

Charging Current: 2A

Operating Power: 10 W

Operating Voltage: DC 3.7 V

Country of origin

Made in: China

Service

Warranty period: 2 years

