

DynaCAD 6.0 hardware and software requirements

Server supported operating system and SQL platforms

Server supported operating system and SQL platforms

Operating systems	SQL Server version
Windows Server 2022*	SQL Server 2019 or 2022
Windows Server 2019	SQL Server 2016, or 2019 or 2022
Windows Server 2016	SQL Server 2016, or 2019 or 2022
Windows 11 Professional, 64-bit*	SQL Server 2019 or 2022
Windows 10 Professional, 64-bit*	SQL Server 2016, or 2019 or 2022

*Windows Server 2022 (Build 20348 or later), Windows 11 (all versions), or Windows 10 (Build 18362/ver. 1903 or later) is required for operation of the QP-Prostate AI DynaCAD option.

The DynaCAD software includes an instance of SQL Express as its database. SQL Standard and Enterprise editions are certified for use with DynaCAD as alternatives but are not included and must be supplied by the customer at the time of installation.

DynaCAD 6.0 servers must meet or exceed the following system configurations

Virtual server		
	Low volume environments where DynaCAD will be processing fewer than 6 cases weekly or connecting less than 4 concurrent users.	High volume environments where DynaCAD will be processing greater than 6 cases weekly or connecting greater than 4 concurrent users – max 10. - or - DynaCAD with QP-Prostate AI option
CPU	8 vCPU minimum	16 vCPU minimum
RAM	32 GB	64 GB minimum
Hard disk	100 GB system + data storage ~300 GB	100 GB system + data storage ~900 GB
Network	1 Gbps bandwidth required*	1 Gbps bandwidth required*
Video card	Onboard integrated card	Onboard integrated card
Display	1280 x 1024 true color	1280 x 1024 true color
VM Environments	DynaCAD can run on verified operating systems on platforms that are certified by Microsoft through the Server Virtualization Validation Program (SVVP)	

*Based on field experience, we have found that the best practice for acceptable system performance is a network latency of <5 ms between server and client.

DynaCAD 6.0 servers must meet or exceed the following system configurations

Physical server hardware

Low volume Environments where DynaCAD will be processing fewer than 6 cases weekly or connecting less than 4 concurrent users.

High volume environments where DynaCAD will be processing greater than 6 cases weekly or connecting greater than 4 concurrent users – max 10.

- or -

DynaCAD with QP-Prostate AI option

CPU	Intel XEON® E5-2600 family, dual processors 4 physical cores per processor with hyperthreading enabled	Intel XEON® E5-2680 family, dual processors 8 physical cores per processor with hyperthreading enabled
Hard disk	100 GB system + data storage ~300 GB	100 GB system + data storage ~900 GB
RAM	32 GB	64 GB
Network	1Gbps bandwidth required*	1Gbps bandwidth required*
Video card	Onboard integrated card	Onboard integrated card
Display	1280 x 1024 true color	1280 x 1024 true color

*Based on field experience, we have found that the best practice for acceptable system performance is a network latency of <5 ms between server and client.

DynaCAD 6.0 clients must meet or exceed the following system configurations

Client supported operating systems

Windows Server 2016, 2019 and 2022

Windows 10 Professional, 64-bit and Windows 11 Professional, 64-bit

Client hardware

CPU	Intel Core® 2 Duo
RAM	8GB
Hard disk	1GB free space
Network	1 Gbps
Video card	Onboard integrated card
Display	1280 x 1024 true color

