




PHILIPS


Phoenix

Phoenix guidewire

Designed to work together
with your rotational
atherectomy system

Phoenix guidewire


Technical specifications



Part number	Type	Diameter	Length	Tip	Style
PG14300LF	Silicone coated nitinol core	.014"	300 cm	Floppy	Light support

Phoenix guidewires are sold in boxes of 5 units each.

Phoenix rotational atherectomy system: a simple and easy to use system



Part number	Catheter size	Introducer size	Working length	Guidewire diameter
P18130K	1,8 mm tracking	5 F (>1,8 mm)	130 cm	0,014"
P18149K	1,8 mm tracking	5 F (>1,8 mm)	149 cm	0,014"
P22130K	2,2mm tracking	6 F (>2,2 mm)	130 cm	0,014"
P22149K	2,2mm tracking	6 F (>2,2 mm)	149 cm	0,014"
PD24127K	2,4 mm deflecting	7 F (>2,4 mm)	127 cm	0,014"