





Product name:	MobileDiagnost wDR	
Identification code(s)	9890-010-89521	
Total weight (in Kg)	555 Kg	
Producer/ Manufacturer	Name company:	Philips Medical Systems DMC GmbH
	Address:	Roentgenstrasse 24
	Zip code:	22335 Hamburg
	Country:	Germany


Recycle Info	Items:	Location
Special attention 	Discharge Power Capacitors before handling the unit	Fig. 4

Fluids / Gases	Items:	Location
	Oil tank / Remove oil before dismantling HV-transformer oil contains no PCBs	Fig. 1,4

Batteries	Type:	Location
	Lead Battery	Fig. 6
	Lithium Battery	Fig. 3 ; 5

Hazardous	Substances:	Location
	Lead	Fig. 1 ; 2

Show locations of materials mentioned on the previous sheet. (Picture information)



Recycling Passport Revision History

Rev	Description	Date
1.0	Initial	9.Dez.2011