





Achieva	Achieva dSTREAM Covers		
Identification code(s)	MAGNET FRONT COVERS dS	9896 030 2564x	
	MAGNET REAR COVERS dS	9896 030 2566x	
	MAGNET SIDE TOP COVERS dS	9896 030 2565x	
Total weight [kg]	± 150 kg		
Manufacturer	Philips Medical Systems Nederland B.V.		
Address	Veenpluis 4-6 Best		
Zip code	5680DA		
Country	The Netherlands		
Electronic info			

Recycling information

Special attention	Items	Location
	Products are made of composite materials consisting of glass, fibers and resins.	
	Microphone printed circuit board in front cover	Figure 1
	4x light sources (LED printed circuit board)	Figure 1

Batteries	Items, type	Location
	Not applicable	

Hazardous substances	Items	Location
	Not applicable	

Fluids and gases	Items	Location
	Not applicable	

Name:	Recycling passport ACHIEVA dSTREAM Covers	12NC:	9896 030 2564x	Version: 00	Doc.Date: 2013-01-30
Author:		Pag.: 1 of 2	Group No: 5790	Section: -	Status: Authorized

Overview of covers:




figure 1

Name:	Recycling passport ACHIEVA dSTREAM Covers	12NC:	9896 030 2564x	Version: 00
Author:				Doc.Date: 2013-01-30
		Pag.: 2 of 2	Group No: 5790	Section: -
				Status: Authorized