

Product name:	dStream Interface 3T		
Identification code(s)	4510-010-1295x		
Total weight (in Kg)	3.8 Kg		
Producer /Manufacturer	Name company:	Philips Healthcare Best	
	Address:	Veenpluis 4-6	
	Zip code:	5680DA Best	
	Country:	The Netherlands	
	Electronic info:	dStream Interface 3T	

Recycle info	Substances:	Location
Special attention 	Items:	
	No special attention needed.	
 Batteries To be removed	Type:	Location
	No batteries present.	
Hazardous  To be removed		Location
	One pcb assembly	Figure 1
	Four so-called RXE modules: 45221330828x or 45980004624x (containing pcb assy)	Figure 2
	Five RF traps containing flex foil pcb's	Figure 3
 Fluids / Gases	Items:	Location
	No fluids and/or gases present.	

Name: WEEE dStream Interface 3T		12NC: DMR122245		Version: 00 Doc.Date: 2010-11-16	
Author:		Pag.: 1 of 3	Group No: 5790	Section: -	Status: see Agile

Locations as mentioned in the Passport (Picture information)




Figure 1, dStream Interface 3T pcb location




Figure 2, dStream Interface 3T RXE locations

Name: WEEE dStream Interface 3T	12NC: DMR122245	Version: 00 Doc.Date: 2010-11-16		
Author:	Pag.: 2 of 3	Group No: 5790	Section: -	Status: see Agile
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2011			Template version 2.	




Figure 3, , dStream Interface 3T RF trap locations

Name: WEEE dStream Interface 3T		12NC: DMR122245		Version: 00 Doc.Date: 2010-11-16	
Author:		Pag.: 3 of 3	Group No: 5790	Section: -	Status: see Agile
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2011				Template version 2.	