

Product name:	dStream Interface 1.5T		
Identification code(s)	4510-010-2487x		
Total weight (in Kg)	3.8 Kg		
Producer /Manufacturer	Name company:	Philips Healthcare Best	
	Address:	Veenpluis 4-6	
	Zip code:	5680DA Best	
	Country:	The Netherlands	
	Electronic info:	dStream Interface 1.5T	

Recycle info	Substances:	Location
Special attention 	Items:	
	No special attention needed.	
 Batteries To be removed	Type:	Location
	No batteries present.	
Hazardous  To be removed		Location
	One pcb assembly	Figure 1
	Four so-called RXE modules: 45221330827x or 45980004622x (containing pcb assy)	Figure 2
	Five RF traps containing flex foil pcb's	Figure 3
 Fluids / Gases	Items:	Location
	No fluids and/or gases present.	

Name: WEEE dStream Interface 1.5T		12NC: DMR122243		Version: 00 Doc.Date: 2010-11-16	
Author:		Pag.: 1 of 3	Group No: 5790	Section: -	Status: see Agile

Locations as mentioned in the Passport (Picture information)




Figure 1, dStream Interface 1.5T pcb location




Figure 2, dStream Interface 1.5T RXE locations

Name: WEEE dStream Interface 1.5T		12NC: DMR122243		Version: 00 Doc.Date: 2010-11-16	
Author:		Pag.: 2 of 3	Group No: 5790	Section: -	Status: see Agile
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2011				Template version 2.	




Figure 3, , dStream Interface 1.5T RF trap locations

Recycling Passport Revision History

Rev	Description	Date
00.1	First Version (Authorised)	21 MAY 2010

Recycling Passport Approval Signatures

Written By	E.Claessens
Approved By	H.d.Zilver

END.

Name: WEEE dStream Interface 1.5T		12NC: DMR122243		Version: 00 Doc.Date: 2010-11-16	
Author:		Pag.: 3 of 3	Group No: 5790	Section: -	Status: see Agile
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2011				Template version 2.	