





<b>Product name</b>	<b>ACC-L 3.0T</b>
Identification code(s)	See list in below tabel
Total weight [kg]	235 ± 20 kg
<b>Manufacturer</b>	<b>Philips Medical Systems Nederland B.V.</b>
Address	<i>Veenpluis 4-6 Best</i>
Zip code	<b>5680DA</b>
Country	<b>The Netherlands</b>
Electronic info	

**Recycling information**

<b>Special Attention</b>		<b>Location</b>
	No Special attention required for the following items	
	1/ ACC2-L Basic Housing (MIU-E only)4598 006 0926x	Figure 1 ACC-L
	3 /Cooling Interface - Assy 4598 005 9555x	Figure 1 ACC-L
	The following items have their own passport and are not part of this passport	<b>Location</b>
	2 /RF AMPLIFIER 3T AN8137 4598 003 2179x	Figure 1 ACC-L

<b>Batteries</b>	<b>Items, type</b>	<b>Location</b>
	Not applicable	

<b>Hazardous substances</b>	<b>Items</b>	<b>Location</b>
	Not Applicable	

<b>Fluids and gases</b>	<b>Items</b>	<b>Location</b>
	Secondary coolant: Distilled or demineralised water with inhibitor and biocide for RF Amplifier circuit Properties: Biocide: 0.01 % ( 3.6ml/40l) of NX1164 Inhibitor: 0.02% (8ml/40l) of AZ8104	Fluid lines

Name:	<b>WEEE RECYCLING PASSPORT &lt;ACC-L&gt;</b>	12NC:	<b>DMR190270</b>	Version: <b>00</b> Doc.Date: <b>2013-11-27</b>
Author:		Page.: <b>1 of 3</b>	Group No: 5790	Section: -
				Status: <b>Authorized</b>

## Locations as mentioned in the Passport

**1- ACC2-L Basic Housing  
(MIU-E only)**

**Mains Inlet unit Extended  
(part of ACC2-L Basic Housing  
cabinet)**

**2 – WC RF AMPLIFIER 3.0T**

**3 – Cooling Interface - Assy**




Figure - 1

Name: <b>WEEE RECYCLING PASSPORT &lt;ACC-L&gt;</b>	12NC: <b>DMR190270</b>	Version: <b>00</b> Doc.Date: <b>2013-11-27</b>
Author:	Page.: <b>2 of 3</b>	Group No: 5790    Section: -
		Status: <b>Authorized</b>

## Identification codes

98960302628x	ACC-L-3T
98960302629x	XR UPGRADE KIT TO 3T
98960302638x	UPGRADE ACC-L CABINET 3T-TX

Name: <b>WEEE RECYCLING PASSPORT &lt;ACC-L&gt;</b>	12NC: <b>DMR190270</b>	Version: <b>00</b> Doc.Date: <b>2013-11-27</b>		
Author:	Page.: <b>3 of 3</b>	Group No: <b>579O</b>	Section: <b>-</b>	Status: <b>Authorized</b>