





Product name	NT-DAC 3T-C-XW6000 / XW6200 / XW6400 / XW6600
Identification code(s)	4522 132 2013x / 4522 132 5079x / 4522 132 5880x / 9896 030 1859x / 9896 030 1821x / 9896 030 2004x / 9896 030 2095x
Total weight [kg]	± 137 kg excluding CDAS (±27kg) and XW6000 /6200 /6400 / 6600 (±17kg)
Manufacturer	Philips Medical Systems Nederland B.V.
Address	Veenpluis 4-6 Best
Zip code	5680DA
Country	The Netherlands
Electronic info	

Recycling information

Special Attention	Items	Location
	NOT part of this passport:	
	1/ CDAS (45221320894x / 45221325996x)	Figure 2
	2/ XW6000/6200(45221322734x /45221324545x)	Figure 3
	XW6200/6400 (45221326520x / 6517x / 6518x /	
	6519x / 7499x) // XW6600 (45221327685x)	
	3/ Stand. Cab. 33HE + 2HE Lade (4522 131 7443 x)	Figure 1


Batteries	Items, type	Location
	Not applicable	

Hazardous substances	Items	Location
	Cable bracket front (Zn yellow passivated 5µm)	Figure 2
	(quick acting) Screw (Zn yellow passivated)	Figure 3
	Cable bracket back (Zn yellow passivated 5µm)	Figure 5
	Cable bracket (Zn yellow passivated 5µm)	Figure 6
	Cable relief bracket (Zn yellow passivated 5µm)	Figure 7

Fluids and gases	Items	Location
	Not applicable	

Name: Recycling Passport <NT-DAC 3T-C-XW6200>	12NC: 4522 132 50790	Version: 02 Doc.Date: 2008-09-04
Author:	Pag.: 1 of 2	Group No: 5790 Section: - Status: Approved

Locations as mentioned in the Passport (Information in pictures)



Name: Recycling Passport <NT-DAC 3T-C-XW6200>		12NC: 4522 132 50790		Version: 02 Doc.Date: 2008-09-04	
Author:	Pag.: 2 of 2	Group No: 5790	Section: -	Status: Approved	Template version 2.
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2008					