



Product name	Solid State RF Amplifier, 15 KW, 64 Mhz
Identification code(s)	4522 132 0324x
Total weight [kg]	81 kg
Manufacturer	Philips Medical Systems Nederland B.V.
Address	Veenpluis 4-6 Best
Zip code	5680DA
Country	The Netherlands
Electronic info	

Recycling information

Special Attention	Items	Location
	None	NA

Batteries	Items, type	Location
	Lithium Ion Button Cell	Figure 2

Hazardous substances	Items	Location
	Printed Circuit Boards Containing Tin / Lead Solder	Figures 1,2,3,4,5

Fluids and gases	Items	Location
	None	NA

Name:	Recycling passport RF Amplifier 1.5T, 53-S30-64	12NC:	4522 132 03240	Version: 00	Doc.Date: 2006-08-07
Author:		Pag.: 1 of 4	Group No: 5790	Section: -	Status: Authorized

Locations as mentioned in the Passport (Information in pictures)





Figure 1
RF Deck – Front View (Power Supply In Rear Not Show)



Note #1
 Printed Circuit Board
 Inside Enclosure



Note #2
 Printed Circuit Board
 Inside Enclosure



Note #3
 Printed Circuit Board
 Inside Enclosure

Name: Recycling passport RF Amplifier 1.5T, 53-S30-64		12NC: 4522 132 03240	Version: 00 Doc.Date: 2006-08-07
Author:	Pag.: 2 of 4	Group No: 5790	Section: - Status: Authorized




Figure 2
Top View




Figure 3
Rear / Top View



Name: Recycling passport RF Amplifier 1.5T, 53-S30-64		12NC: 4522 132 03240		Version: 00 Doc.Date: 2006-08-07	
Author:		Pag.: 3 of 4	Group No: 5790	Section: -	Status: Authorized




Figure 4
Power Supply – Left Side
 (Removed From Enclosure For Clarity)



Figure 5
Power Supply – Right Side
 (Removed From Enclosure For Clarity)

Name: Recycling passport RF Amplifier 1.5T, 53-S30-64		12NC: 4522 132 03240		Version: 00 Doc.Date: 2006-08-07	
Author:		Pag.: 4 of 4	Group No: 5790	Section: -	Status: Authorized