





<b>Product name</b>	<b>WA Busbars</b>
Identification code(s)	Single mode 989603023064 Dual mode 989603023154
Total weight [kg]	14 kg
<b>Manufacturer</b>	<b>Philips Medical Systems Nederland B.V.</b>
Address	Veenpluis 4-6 Best
Zip code	5680DA
Country	The Netherlands
Electronic info	

## Recycling information

Special Attention	Items	Location
	Copper bars are wrapped and injected with composite materials consisting of glass, fibers and epoxy resins.	

Batteries	Items, type	Location
	N.A.	

Hazardous substances	Items	Location
	N.A.	

Fluids and gases	Items	Location
	N.A.	

Name: <b>Recycling Passport WA BusBars</b>		12NC: <b>9896030230641</b>		Version: <b>01</b> Doc.Date: <b>2011-09-26</b>	
Author: <b>Menno Mathlener</b>		Pag.: <b>1 of 1</b>	Group No: <b>5790</b>	Section: <b>-</b>	Status: <b>Authorized</b>