





Product name	Gradient Coil: WA30S
Identification code(s)	4522 133 0576.x
Total weight [kg]	< 850 kg (typical 825 kg)
Manufacturer	Tesla Engineering Limited
Address	Water Lane, Storrington, Sussex
Zip code	RH20 3EA
Country	UK
Electronic info	

Recycling information

Special Attention	Items	Location
	Coil contains cooling water + additives, see below.	Cooling circuit, see Figure 1

Batteries	Items, type	Location
	No batteries	

Hazardous substances	Items	Location
	PBB in FR4 insulation sheets	Between copper coil layers and copper coil-to-coil connection strips (see Figure 2)

Fluids and gases	Items	Location
	Additives to cooling water:	Cooling circuit, see Figure 1
	Betz Dearborn Inhibitor AZ8104 ; 0.2mL/L	
	Betz Dearborn Biocide Spectrus NX1164, 0.09mL/L or Betz Dearborn Biocide Spectrus NX1106, 0.09mL/L	

Name: Recycling passport WA30S Gradient Coil		12NC: 4522 133 05760		Version: R0 V02	
Author:		Pag.: 1 of 2	Group No: 5790	Section: 00	Doc.Date: 2011-06-24
			Status: Authorized		

Locations as mentioned in the Passport (Information in pictures)

Figure 1.




Figure 2.

Name: Recycling passport WA30S Gradient Coil			12NC: 4522 133 05760		Version: R0 V02 Doc.Date: 2011-06-24	
Author:		Pag.: 2 of 2	Group No: 5790	Section: 00	Status: Authorized	