

Product name:	Anterior Coil 3T		
Identification code(s)	4598-001-7120x		
Total weight (in Kg)	4.5 Kg		
Producer /Manufacturer	Name company:	Philips Healthcare Best	
	Address:	Veenpluis 4-6	
	Zip code:	5680DA Best	
	Country:	The Netherlands	
	Electronic info:	Anterior Coil 3T	

Recycle info	Substances:	Location
Special attention 	Items:	
	No special attention needed.	
 Batteries To be removed	Type:	Location
	No batteries present.	
Hazardous  To be removed		Location
	Sixteen pcb assemblies and one Flex-foil pcb assembly	Figure 1
	Four so-called RXE modules: 45980004624x (containing pcb assy)	Figure 2
 Fluids / Gases	Items:	Location
	No fluids and/or gases present.	

Name:	WEEE Anterior Coil 3T	12NC:	DMR147880	Version: 01
Author:		Group No: 5790	Section: -	Doc.Date: 2012-05-14
	Pag.: 1 of 3			Status: see Agile

Locations as mentioned in the Passport (Picture information)




Figure 1, Anterior Coil 3T pcb locations

Name: WEEE Anterior Coil 3T		12NC: DMR147880		Version: 01 Doc.Date: 2012-05-14	
Author:		Pag.: 2 of 3	Group No: 5790	Section: -	Status: see Agile
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2012				Template version 2.	




Figure 2, Anterior Coil 3T RXE locations

Name: WEEE Anterior Coil 3T		12NC: DMR147880		Version: 01 Doc.Date: 2012-05-14	
Author:		Pag.: 3 of 3	Group No: 5790	Section: -	Status: see Agile
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2012				Template version 2.	