






PHILIPS HEALTHCARE	Product Recycling Passport	Page 1 of 2
-----------------------	-----------------------------------	-------------

Product name	Module Charger (REF 989803191031)	
Identification code(s)	REF 989803191031	
Total weight (in Kg)	0.38 kg (without packaging)	
Producer/ Manufacturer	Name company:	Philips HealthTech
	Address:	3000 Minuteman Rd. Andover, MA
	Zip code:	01810
	Country:	U.S.A.

Recycle Info	Items:	Location
Special attention 	Packaging (Corrugated Board) (total weight = 0.42kg)	Packaging
	Plastic (PC/ABS blend, total weight = 0.29 kg)	Outer Enclosure
	PCB (total weight = 0.078 kg)	Internal PCB
	Copper (total weight = 0.0018 kg)	Cables
 Fluids / Gases	Items:	Location
	N/A	N/A
Batteries  To be Removed	Type:	Location
	N/A	N/A
Hazardous  To be Removed	Substances:	Location
	N/A	N/A

Show locations of materials mentioned on the previous sheet.
 (Please use an illustration/drawing of the system, photos are discouraged)



Recycling Passport Revision History

Rev	Description	Date
A	Initial Release	20-Jun-2016

Recycling Passport Approval Signatures

(Note: Do not include the section below in the published version to protect privacy)