





PHILIPS HEALTHCARE	Product Recycling Passport	<b>PHILIPS</b>
--------------------	----------------------------	----------------

<b>Product name:</b>	MRC200 0407 ROT-GS 1004   MRC200 0508 ROT-GS 1003 MRC200+ 0407 ROT-GS 1004   MRC200+ 0508 ROT-GS 1003
<b>Identification code(s):</b>	9890-000-86501   9890-000-86521 9890-000-85103   9890-000-85143 9890-000-86741   9890-000-86751
<b>Total weight (in kg):</b>	75,5 kg   67 kg

<b>Producer / Manufacturer</b>	<b>Name company:</b>	Philips Medical Systems DMC GmbH
	<b>Address:</b>	Röntgenstrasse 22
	<b>Zip code:</b>	22335 Hamburg
	<b>Country:</b>	Germany

Recycle Info		
<b>General Hazards</b> <i>(Desinfection)</i>	<b>Items</b>	<b>Location</b>
	<i>For dismantling activities, Treatment Facilities, must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account.</i>	
<b>Special attention</b>  	<b>Items</b>	<b>Location</b>
 <b>Fluids/ Gases</b>	<b>Items</b>	<b>Location</b>
	Oil: 5,4 kg for 9890-000-86501 & 9890-000-85103	
	Oil: 5,4 kg for 9890-000-86741	
	Oil: 5 kg for 9890-000-86521 & 9890-000-85143	
	Oil: 5 kg for 9890-000-86751	
 <b>Batteries</b> <i>to be removed</i>	<b>Type</b>	<b>Location</b>
 <b>Hazardous</b> <i>to be removed</i>	<b>Substances</b>	<b>Location</b>
	Lead as x-ray shielding inside housing: 4,8 kg	Figure 1

Printed copies are uncontrolled duplicates

GTC- Template Recycling Passport XMT-GTC-03-0196 Vers. 1.0

**Location of materials**

Lead as x-ray shield




Figure 1

**Revision History**

Rev	Description	Date
1.0	First signed version	03. Nov 2015

Printed copies are uncontrolled duplicates

GTC- Template Recycling Passport XMT-GTC-03-0196 Vers. 1.0