





Product name:	QBC Panorama 1T	
Identification code(s)	Product number	4522 139 00190
	Consists of:	
	Discs (2 pcs)	4522 132 3245X
	Cover disc (plastic) (2pcs)	4522 132 4271X
	Box (brass)	4522 137 0030X
	Box (brass)	4522 137 0029X
	Box cable assembly	4522 137 0034X
	Rigid cable sets	4522 132 3251X
		4522 132 3250X
		4522 132 3252X
	Cooling assembly with resistor	4522 132 3683X
Total weight (in Kg)	48 Kg	

Producer /Manufacturer	Name company:	Philips Medical Systems Best
	Address:	Veenpluis 4-6
	Zip code:	5680DA Best
	Country:	The Netherlands
	Electronic info:	RF_Panorama1TQBC_452213900190_579_v01.doc

Recycle info	Substances:	Location
Special attention 	Items:	
	No special attention needed.	

Batteries	Type:	Location
 To be removed	No batteries present.	

Hazardous		Location
 To be removed	Cooling assembly with resistor 4522 132 3683X contains Beryllium-oxyde	Figure 1
	Printed Circuit Board Assemblies are present in:	
	Disc (6 Printed Circuit Board Assemblies) 4522 132 3245X	Figure 2
	Box (brass) (2 Printed Circuit Board Assemblies) 4522 137 0030X	Figure 3
	Box (brass) (1 Printed Circuit Board Assembly) 4522 137 0029X	Figure 4
	Box cable assembly (2 Printed Circuit Board Assemblies) 4522 137 0034X	Figure 5

	Items:	Location
	No fluids and/or gasses present.	

Name:	WEEE passport QBC Panorama 1T	12NC:	4522 139 00190	Version: 01
Fluids / Gases				Doc.Date: 14-08-2006
Author: M. Vrinten		Pag.: 1 of 4	Group No: 579	Section: 0
				Status: authorized

Locations as mentioned in the Passport (Pictures information)

Figure 1, Cooling assembly with resistor (4522 132 3683X) containing beryllium-oxyde




Figure 2, Discs (4522 132 3245X), each containing six Printed Circuit Board Assemblies





Figure 3, Brass box (4522 137 0030X) containing two Printed Circuit Board Assemblies



Name: WEEE passport QBC Panorama 1T		12NC: 4522 139 00190	Version: 01 Doc.Date: 14-08-2006	
Author: M. Vrinten	Pag.: 2 of 4	Group No: 579	Section: 0	Status: authorized

Figure 4, Brass box (4522 137 0030X) containing one Printed Circuit BoardA





Figure 5, Box cable assembly (4522 137 0034X) containing two Printed Circuit Board's



Name:	WEEE passport QBC Panorama 1T	12NC:	4522 139 00190	Version: 01
Author: M. Vrinten		Pag.: 3 of 4	Group No: 579	Section: 0
				Doc.Date: 14-08-2006
				Status: authorized

Rev	Description	Date
1	First Passport	August 9, 2006

Recycling Passport Approval Signatures

Written By

M. Vrinten

Reviewed By (optional)

A. Snelders

Reviewed By (optional)

J. Jonkergouw

Approved By

GJ Plattel

Name: WEEE passport QBC Panorama 1T		12NC: 4522 139 00190	Version: 01 Doc.Date: 14-08-2006
Author: M. Vrinten	Pag.: 4 of 4	Group No: 579	Section: 0 Status: authorized