




Product name	Invivo Sense Knee coils
Identification code(s)	4522 132 31070 8-SKC 1.5T 4535 670 61680 8-SKC 3.0T 9896 030 13380 6-SKC 3.0T
Total weight [kg]	8.5 kg
Manufacturer	Invivo, corporation. USA
Address	3545 SW 47th Avenue
Zip code	Gainesville, FL 32608
Country	USA
Electronic info	RFCoils_Invivo_SenseKneeCoils.doc

Recycling information

Special Attention	Items	Location
	None	Figure

Batteries	Items, type	Location
	None	Figure

Hazardous substances	Items	Location
	Several printed circuit boards inside the coil	Figure 1
	Beryllium copper contact pins on connectors	Figure 1

Fluids and gases	Items	Location
	6-ch 3T coil: Plastic bottle with copper sulfate	
	Bottle is delivered with the coil separately	

Name Recycling passport Invivo SKC	12NC: 4522 132 31070	Version: 01 Doc.Date: 2007-05-11
Author: Paul Ebben	Page 1 of 3	Group No: 5790 Section: -
		Status: Authorized

Locations as mentioned in the Passport (Information in pictures)

Place sketch (es) or photo(s) here, and indicate issues using icons:



PWB's are inside the coil

Beryllium contacts inside the coil + in the connector.



Figure 1

Name Recycling passport Invivo SKC	12NC: 4522 132 31070	Version: 01 Doc.Date: 2007-05-11
Author: Paul Ebben	Page 2 of 3	Group No: 5790 Section: -
		Status: Authorized

Recycling Passport Revision History

Revision	Description	Date
00	Initial version	2006-08-04

Recycling Passport Approval Signatures

<p>Written By Paul Ebben</p>	<p>Reviewed By (optional) Geert-Jan Platel Eugene Fatica</p>
<p>Approved By Paul Ebben</p>	<p>Approved By (optional)</p>

End of document.

Name Recycling passport Invivo SKC	12NC: 4522 132 31070	Version: 01 Doc.Date: 2007-05-11
Author: Paul Ebben	Page 3 of 3	Group No: 5790 Section: -
		Status: Authorized