





Product name	Invivo Knee/Foot coil
Identification code(s)	4522 132 14020 Knee Foot 1.5T
Total weight [kg]	4 kg
Manufacturer	Invivo, corporation. USA
Address	3545 SW 47th Avenue
Zip code	Gainesville, FL 32608
Country	USA
Electronic info	RFCoils_Invivo_KneeFootCoil.doc

Recycling information

Special Attention 	Items	Location
	None	Figure

Batteries 	Items, type	Location
	None	Figure

Hazardous substances 	Items	Location
	Beryllium copper contact pins on connectors	Figure 1
	Several printed circuit boards inside the coil	Figure 1

Fluids and gases 	Items	Location
	1.5T coil → plastic bottle with copper sulfate	See next page
	0.008moles/liter CuSO ₄ and 0.034moles/liter	

Name: Recycling passport Invivo Knee/Foot coil		12NC: 4522 150 14020		Version: 00 Doc.Date: 2006-08-04	
Author: Paul Ebben		Page 1 of 3	Group No: 5790	Section: -	Status: Authorized

Locations as mentioned in the Passport (Information in pictures)




Figure 1

Name: Recycling passport Invivo Knee/Foot coil	12NC: 4522 150 14020	Version: 00 Doc.Date: 2006-08-04
Author: Paul Ebben	Page 2 of 3	Group No: 5790 Section: -
		Status: Authorized

Recycling Passport Revision History

Revision	Description	Date
00	Initial version	2006-08-04

Recycling Passport Approval Signatures

<p>Written By Paul Ebben</p>	<p>Reviewed By (optional) Geert Jan Platel Bert Jansen</p>
<p>Approved By Paul Ebben</p>	<p>Approved By</p>

End of document.

Name: Recycling passport Invivo Knee/Foot coil	12NC: 4522 150 14020	Version: 00 Doc.Date: 2006-08-04
Author: Paul Ebben	Page 3 of 3	Group No: 5790 Section: -
		Status: Authorized