





Product name	GCI
Identification code(s)	4522-131-9114x GCI Subassy 4522-117-9275x GCI1 Mainboard 4522-117-9529x GCI2 Mainboard 4522-117-9280x GCI1 IO board 4522-117-9496x GCI2_IO_board
Total weight [kg]	5 kg
Manufacturer	Philips Medical Systems Nederland B.V.
Address	Veenpluis 4-6 Best
Zip code	5680DA
Country	The Netherlands
Electronic info	GradientInterfaceGCI_452213191140_579_v01.doc

Recycling information

Special Attention	Items	Location
	PVC wiring & cabling	Figure 1

Batteries	Items, type	Location
	---	

Hazardous substances	Items	Location
	Printed circuit boards utilize lead	Internal to assembly
	Be-Cu contact springs on front panel	Front panel, Figure 2

Fluids and gases	Items	Location
	---	Figure ..

Name: Recycling passport GCI		12NC: 4522 131 9114 x		Version: 01 Doc.Date: 2006-09-11	
Author: Martin Hollander		Page 1 of 3	Group No: 5790	Section: -	Status: Authorized

Locations as mentioned in the Passport (Information in pictures)

Place sketch(es) or photo(s) here, and indicate issues using icons:




Figure 1



Figure 2

Name: Recycling passport GCI		12NC: 4522 131 9114 x	Version: 01 Doc.Date: 2006-09-11
Author: Martin Hollander	Page 2 of 3	Group No: 5790	Section: -
			Status: Authorized

Recycling Passport Revision History

Revision	Description	Date
1	Initial version	2006-Sept-11

Recycling Passport Approval Signatures

Written By Martin Hollander	Reviewed By (optional)
Approved By Nico Swennen	Approved By

End of document.

Name: Recycling passport GCI		12NC: 4522 131 9114 x		Version: 01 Doc.Date: 2006-09-11	
Author: Martin Hollander		Page 3 of 3	Group No: 5790	Section: -	Status: Authorized