

Product name	COMFORT ZONE
Identification code(s)	9896 030 0744x
Total weight [kg]	10.6 kg (including packaging)
Manufacturer	Philips Medical Systems Nederland B.V.
Address	Veenpluis 4-6 Best
Zip code	5680DA
Country	The Netherlands
Electronic info	ComfortZone_989603007440_579_v01.doc

Recycling information

Special Attention 	Items	Location
	The 3 mirrors are made of glass.	Figure 3

Batteries 	Items, type	Location
	No	

Hazardous substances 	Items	Location
		Figure 3,4,5

Fluids and gases 	Items	Location
	No	

Name: Recycling passport COMFORT ZONE				12NC: 9896 030 0744 x		Version: 00 Doc.Date: 2006-07-19	
Author: D.Sinnema		Page 1 of 4		Group No: 5790	Section: -	Status: Released	

Locations as mentioned in the Passport (Information in pictures)

Figure 1: complete Comfort Zone option



Materials used:

Comfort Zone stand including small camera (incl PCB), 3 mirrors, cable

The housings are made of Polycarbonate

CZ-interface (PCB)

Cable clip (copper)

Filterbox including PCB

Videofilter including PCB

Cables

Name:	Recycling passport COMFORT ZONE			12NC:	9896 030 0744 x		Version: 00
Author: D.Sinnema			Page 2 of 4	Group No: 5790	Section: -		Doc.Date: 2006-07-19
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2006				Status: Released			

Figure 2: Comfort Zone stand



Figure 3: mirrors highlighted



Figure 4: video switch (including PCB inside)



Name: Recycling passport COMFORT ZONE				12NC: 9896 030 0744 x		Version: 00 Doc.Date: 2006-07-19	
Author: D.Sinnema			Page 3 of 4	Group No: 5790	Section: -	Status: Released	

© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2006

Figure 5: Comfort Zone PCB and filterbox (including PCB inside)



Recycling Passport Revision History

Revision	Description	Date
1	Initial version	19-07-2006

Recycling Passport Approval Signatures

Written By D.Sinnema	Reviewed By (optional)
Approved By D.Sinnema	Approved By

End of document.

Name: Recycling passport COMFORT ZONE	12NC: 9896 030 0744 x	Version: 00 Doc.Date: 2006-07-19
Author: D.Sinnema	Page 4 of 4	Group No: 5790 Section: - Status: Released