

Recycling Passport

Product name:	BuckyDiagnost TH2/TF		
Identification code	9890-010-06503		
Order No.			
Device Type	B2B (Business to business)		
Producer /Manufacturer	Company name:	Philips Medical Systems DMC GmbH	
	Address:	Roentgenstrasse 24	
	Zip code:	22335 Hamburg	
	Country:	Germany	
	Web info:	See Philips Intranet	
Total weight (in Kg) thereof			183 Kg – 289 Kg

Version: TH2 (750 table width; height adjustable)			
Klasse	Merkmalbezeichnung	Gesamtgewicht	Gewichtseinheit
steel, iron	iron, low alloy (<5%)	171,852	KG
	iron, high alloy (>5%)	17,336	KG
steel, iron		189,188	KG
nonferrous metals and alloys	aluminum, -alloy	3,634	KG
	copper, -alloy	0,001	KG
	zinc, -alloy	0,001	KG
nonferrous metals and alloys		3,636	KG
glass / ceramics	ferrite	0,025	KG
glass / ceramics		0,025	KG
plastics / organic substances	elastomer	8,169	KG
	thermoplastic	2,43	KG
	lubricating grease	0,007	KG
	thermoset	0,002	KG
plastics / organic substances		10,608	KG
batteries	batteries, lithium	0,012	KG
batteries		0,012	KG
standard parts	bearing plates	43	KG
	cables	11,469	KG
	motors, pumps	10,299	KG
	rasters	2,23	KG
	power supplies / transformers	1,98	KG
	wheels / castors	1,3	KG
	brakes / clutches	0,83	KG
	other electronic powered devic	0,804	KG
	printed circuit boards	0,419	KG
	mounting parts, attaching part	0,084	KG
standard parts		72,415	KG
relevant materials	confirmed by supplier		
relevant materials			
prohibited materials	confirmed by supplier		
prohibited materials			
		275,884	KG

Printed copies are uncontrolled documents

Recycling Passport

Version:TH2 (850 table width; height adjustable)

Klasse	Merkmalbezeichnung	Gesamtgewicht	Gewichtseinheit
steel, iron	iron, low alloy (<5%)	173,952	KG
	iron, high alloy (>5%)	17,336	KG
steel, iron		191,288	KG
nonferrous metals and alloys	aluminum, -alloy	3,634	KG
	copper, -alloy	0,001	KG
	zinc, -alloy	0,001	KG
nonferrous metals and alloys		3,636	KG
glass / ceramics	ferrite	0,025	KG
glass / ceramics		0,025	KG
plastics / organic substances	elastomer	8,169	KG
	thermoplastic	2,43	KG
	lubricating grease	0,007	KG
	thermoset	0,002	KG
plastics / organic substances		10,608	KG
batteries	batteries, lithium	0,012	KG
batteries		0,012	KG
standard parts	patient tables	50,5	KG
	cables	11,469	KG
	motors, pumps	10,299	KG
	rasters	2,23	KG
	power supplies / transformers	1,98	KG
	wheels / castors	1,3	KG
	brakes / clutches	0,83	KG
	other electronic powered devic	0,804	KG
	printed circuit boards	0,419	KG
mounting parts, attaching part	0,084	KG	
standard parts		79,915	KG
relevant materials	confirmed by supplier		
relevant materials			
prohibited materials	confirmed by supplier		
prohibited materials			
		285,484	KG

Printed copies are uncontrolled documents

Recycling Passport

Version:TF (750 table width; height fixed)





Klasse	Merkmalbezeichnung	Gesamtgewicht	Gewichtseinheit
steel, iron	iron, low alloy (<5%)	95,352	KG
	iron, high alloy (>5%)	17,336	KG
steel, iron		112,688	KG
nonferrous metals and alloys	aluminum, -alloy	3,684	KG
	copper, -alloy	0,005	KG
	zinc, -alloy	0	KG
nonferrous metals and alloys		3,689	KG
glass / ceramics	ferrite	0,025	KG
glass / ceramics		0,025	KG
plastics / organic substances	thermoplastic	2,789	KG
	elastomer	2,088	KG
	lubricating grease	0,006	KG
	thermoset	0,002	KG
plastics / organic substances		4,885	KG
batteries	batteries, lithium	0,012	KG
batteries		0,012	KG
standard parts	bearing plates	43	KG
	cables	11,469	KG
	rasters	2,23	KG
	power supplies / transformers	1,98	KG
	wheels / castors	0,98	KG
	brakes / clutches	0,83	KG
	other electronic powered devic	0,804	KG
	printed circuit boards	0,708	KG
	motors, pumps	0,099	KG
	mounting parts, attaching part	0,084	KG
standard parts		62,184	KG
relevant materials	confirmed by supplier		
relevant materials			
prohibited materials	confirmed by supplier		
prohibited materials			
		183,483	KG

Printed copies are uncontrolled documents

Option:Tomo (motor)

Klasse	Merkmalbezeichnung	Gesamtgewicht	Gewichtseinheit
steel, iron	iron, low alloy (<5%)	0,783	KG
steel, iron		0,783	KG
nonferrous metals and alloys	aluminum, -alloy	0,15	KG
nonferrous metals and alloys		0,15	KG
plastics / organic substances	thermoplastic	0,057	KG
plastics / organic substances		0,057	KG
standard parts	motors, pumps	1,75	KG
	cables	0,703	KG
	mounting parts, attaching part	0,05	KG
	other electronic powered devic	0,001	KG
standard parts		2,504	KG
		3,494	KG

Recycling Passport

Recycle info		
General Hazards (Desinfection)	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
Hazardous Substances  To be removed	Substances:	Location
Batteries  To be removed	Type:	Location
	batteries, lithium	Fig. 1
Special attention 	Item	Location
Fluids / Gases 	Item	Location
	HV-transformer oil contains no PCBs	

Printed copies are uncontrolled documents

Recycling Passport

Locations as mentioned in the Passport (Pictures information)

Version: TH2 (750 or 850 table width; height adjustable)



1



if sensing

Version: TF (750 table width; height fixed)



Printed copies are
uncontrolled documents