





Recycling Passport

Product name:	ROT 501 / ROT-GS 501 / ROT-GS 505 / ROT 507 SRM xxxx		
Identification code	9874-008-64132 9890-000-03621 9890-000-63912 9890-000-63671 9890-000-63721 9874-008-46122 9890-000-03362	9890-000-03822 9890-000-85182	9890-000-63941 9890-000-03442 9890-000-63931 9890-000-03562
Order No.			
Device Type	B2B (Business to business)		
Producer /Manufacturer	Company name:	Philips Medical Systems DMC GmbH	
	Address:	Roentgenstrasse 24	
	Zip code:	22335 Hamburg	
	Country:	Germany	
	Web info:	http://staging.medical.philips.com/main/company/sustainability/recycling/	
Total weight (in Kg) thereof			27,5

steel, iron	iron, low alloy (<5%)	3,1KG
	iron, high alloy (>5%)	0,7KG
steel, iron		3,8KG
nonferrous metals and alloys	aluminium and -alloys	6,8KG
	copper and -alloys	2,7KG
	molybdenum and -alloys	0,9KG
nonferrous metals and alloys		10,4KG
glass / ceramics	ferrite	0,3KG
	glass	0,1KG
glass / ceramics		0,4KG
plastics / organic substances	oil	4,2KG
	thermoset	2KG
	elastomer	0,1KG
plastics / organic substances		6,3KG
relevant materials	lead and -alloys	5,5KG
	tungsten and -alloys	0,1KG
	beryllium and -alloys	0,001KG
relevant materials		5,601KG

Printed copies are uncontrolled documents

Recycling Passport

Recycle info		
General Hazards (Desinfection)	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
Hazardous Substances  To be removed	Substances:	Location
	Lead as x-ray shielding inside housing: 5,5 kg	Figure 1
	Beryllium as x-ray window in vacuum envelope: 1 g	Figure 1
Batteries  To be removed	Type:	Location
Special attention 	Item	Location
Fluids / Gases 	Item	Location
	Oil: 4,2 kg	
	HV-transformer oil contains no PCBs	

Printed copies are uncontrolled documents

Recycling Passport

Locations as mentioned in the Passport (Pictures information)




Figure 1

Printed copies are
uncontrolled documents