





## Recycling Passport

<b>Product name:</b>	Velara R/F 2T	Velara R/F 2T	Velara R/F Rel.2
<b>Identification code</b>	9890-000-62062	9890-000-62063	9890-000-70151
<b>Product name:</b>	Velara R/F 2T Rel.3	Velara R/F 2T Rel.3	
<b>Identification code</b>	9890-000-70291	9890-000-70294	
<b>Product name:</b>	Velara R/F 2T R.3.11.1		
<b>Identification code</b>	9890-000-70295		
<b>Order No.</b>			
<b>Device Type</b>	B2B (Business to business)		
<b>Producer /Manufacturer</b>	Company name:	Philips Medical Systems DMC GmbH	
	Address:	Roentgenstrasse 24	
	Zip code:	22335 Hamburg	
	Country:	Germany	
	Web info:	<a href="http://www.medical.philips.com/main/company/sustainability/recycling/xray_generators/">http://www.medical.philips.com/main/company/sustainability/recycling/xray_generators/</a>	
<b>Total weight (in Kg)</b>			285
<b>thereof</b>			

steel, iron	iron, low alloy (<5%)	130,785	KG
	iron, high alloy (>5%)	2,334	KG
<b>steel, iron</b>		<b>133,119</b>	<b>KG</b>
nonferrous metals and alloys	aluminum, -alloy	12,274	KG
	copper, -alloy	5,005	KG
	zinc, -alloy	1,576	KG
<b>nonferrous metals and alloys</b>		<b>18,855</b>	<b>KG</b>
glass / ceramics	glass	0,288	KG
<b>glass / ceramics</b>		<b>0,288</b>	<b>KG</b>
plastics / organic substances	oil	48,06	KG
	thermoplastic	5,952	KG
	thermoset	4,245	KG
	elastomer	0,232	KG
<b>plastics / organic substances</b>		<b>58,489</b>	<b>KG</b>
standard parts	other electronic powered devices	45,92	KG
	printed circuit boards	11,584	KG
	cables	6,733	KG
	power supplies / transformers	6,184	KG
	mounting parts, attaching part	2,03	KG
	fans	1,804	KG
	motors, pumps	0,52	KG
<b>standard parts</b>		<b>74,775</b>	<b>KG</b>
		<b>285,526</b>	<b>KG</b>

# Recycling Passport

<b>Recycle info</b>		
<b>General Hazards (Desinfection)</b>	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
<b>Hazardous Substances</b>    <b>To be removed</b>	<b>Substances:</b>	<b>Location</b>
<b>Batteries</b>    <b>To be removed</b>	<b>Type:</b>	<b>Location</b>
	Lithium chrome cell, 3V	Figure 1 (1)
	Lithium chrome cell, 3V	Figure 1 (2)
<b>Special attention</b>  	<b>Item</b>	<b>Location</b>
	Pb is present in the soldering process of PCBs	Figure 1
<b>Fluids / Gases</b>  	<b>Item</b>	<b>Location</b>
	Transformer Oil, Shell Diala	Figure 1 (3)
	HV-transformer oil contains no PCBs	

# Recycling Passport




Figure 1