





Recycling Passport

Product name:	Cooling unit CU 3101		
Identification code	9890 000 85652		
Order No.			
Device Type	B2B (Business to business)		
Producer /Manufacturer	Company name:	Philips Medical Systems DMC GmbH	
	Address:	Roentgenstrasse 24	
	Zip code:	22335 Hamburg	
	Country:	Germany	
	Web info:	http://staging.medical.philips.com/main/company/sustainability/recycling/	
Total weight (in Kg) thereof			37,0

steel, iron	iron, low alloy (<5%)	12	KG
	iron, high alloy (>5%)	3,7	KG
steel, iron		15,7	KG
nonferrous metals and alloys	aluminum and -alloys	12	KG
	copper and -alloys	2,8	KG
nonferrous metals and alloys		14,8	KG
plastics / organic substances	oil	5	KG
	polyamide	1,3	KG
	elastomer	0,2	KG
plastics / organic substances		6,5	KG

Printed copies are uncontrolled documents

Recycling Passport

Recycle info		
General Hazards (Desinfection)	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
Hazardous Substances  To be removed	Substances:	Location
Batteries  To be removed	Type:	Location
Special attention 	Item	Location
Fluids / Gases 	Item	Location
	Oil: 5,0 kg	
	HV-transformer oil contains no PCBs	

Printed copies are uncontrolled documents

Recycling Passport

Locations as mentioned in the Passport (Pictures information)

Printed copies are
uncontrolled documents



Figure 1