





Product name	NTDAC 1.5T XW6000 / (NTDAC 1.0T / NTDAC 0.5T)
Identification code(s)	45221324445x / 45221322018x / 45221321957x / 45221320015x / 45221321887x / 45221320014x / 45221321885x / 45221321955x / 45221322016x / 45221321883x / 45221320013x / 45221321953x
Total weight [kg]	137 kg Excluding BDAS (± 48 kg) and W6000 / XW6000 (± 16 kg)
Manufacturer	Philips Medical Systems Nederland B.V.
Address	Veenpluis 4-6, Best
Zip code	5684 PC
Country	The Netherlands
Electronic info	System_NTDAC1.5TforBDAS_452213244450_579_v01

Recycling information

Special Attention	Items NOT part of this passport:	Location
	1/ BDAS (45221308900x)	Figure 1
	2/ W6000 (45221321110x) / W6000 (45221321951x)	Figure 1
	XW6000 (45221321973x) / XW6000 (45221322733x)	Figure 1
	3/ STAND.CAB. 33E + 2HE LADE (45221317443x)	Figure 1

Batteries	Items, type	Location
	Not applicable	


Hazardous substances	Items	Location
	(quick acting) Screw (Zn yellow passivated 5µm)	Figure 3
	Separating plate (Zn yellow passivated 5µm)(if applicable)	Figure 4
	Cable relief bracket (Zn yellow passivated 5µm)	Figure 5

Fluids and gases	Items	Location
	Not applicable	

Name: Recycling Passport NTDAC 1.5T XW6000		12NC: 4522 132 44450	Version: 01 Doc.Date: 2007-07-12
Author: G. Geboers	Pag.: 1 of 3	Group No: 5790	Section: - Status: Authorized

Locations as mentioned in the Passport (Information in pictures)

<Place sketch(es) or photo(s) here, and indicate issues using icons >



BDAS

W 6000/ XW6000

Stand. Cabinet



Figure 3



Figure 1: Frontview
NT-DAC cabinet




Figure 4



Figure 2: Backview
NT-DAC cabinet



Figure 5



Name: Recycling Passport NTDAC 1.5T XW6000		12NC: 4522 132 44450		Version: 01 Doc.Date: 2007-07-12	
Author: G. Geboers		Pag.: 2 of 3	Group No: 5790	Section: -	Status: Authorized

Recycling Passport Revision History

Revision	Description	Date
1	Initial version	2007-07-12

Recycling Passport Approval Signatures

Written By G. Geboers	Reviewed By (optional)
Approved By Geert-Jan Plattèl	Approved By

<END OF DOCUMENT>

Name: Recycling Passport NTDAC 1.5T XW6000	12NC: 4522 132 44450	Version: 01 Doc.Date: 2007-07-12		
Author: G. Geboers	Pag.: 3 of 3	Group No: 5790	Section: -	Status: Authorized