





Product name:	7T Multinuclear Partial Volume Coil (C13/1H)	
Identification code	4535 674 5756*	
Total weight (in Kg)	Approximately 2.8 Kg	
Producer /Manufacturer	Name company:	Philips Medical Systems (Cleveland), Inc
	Address:	595 Miner Rd., Cleveland, Ohio
	Zip code:	44143
	Country:	USA
	Electronic info:	http://www.medical.philips.com/main/company/sustainability/recycling/mri/

Recycle info	Substances	Location
Hazardous  To be removed	Lead (Pb) (Solder - PC boards)	Figure 1
 Batteries To be removed	Type	Location
Special attention 	Item	Location
 Fluids / Gases	Item	Location

* Denotes revision level wildcard designation. (Numbers range from 1-9).

Locations as mentioned in the Passport (Pictures information)




Figure 1 7T Multinuclear 13C/1H partial volume coil.